

Raxwell



RX Clear Vinyl Gloves



Powder-Free
 Latex-Free
 100 Gloves/Box
 Chemical & Oil Resistance
 Multi-Purpose
 Top Quality

Standards

(US) FDA 510K / ASTM D6319 / ASTM D6978-05(2019) / ASTM F2671-15

(EU) CE / EN 455 part 1, 2, 3&5 / EN 374 part 1, 2, 3&4 / EN 16523-1:2015 / EN ISO 21420-2020

(CN) GB4806.11-2016

Particular	Specifications	
Palm width (mm)	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X Large	125 ± 5
Length (mm)	235 ± 5	
Thickness (Single Wall)	3 mil	
Tensile Strength	Before aging (MPa)	min. 13
	After aging (MPa)	min. 12
Elongation	Before aging (%)	400
	After aging (%)	380

Part No.	Size	LBS /Carton	Cartons/Pallet
RX60001	S	8.91	84
RX60002	M	9.83	84
RX60003	L	10.96	84
RX60004	XL	11.95	84
Packing	100Pcs/Color Box By Weight, 10Boxes/CTN		
Box Dimension	8.27x4.53x2.52in		
Carton Dimension	12.99x9.45x8.66in		
Contain Loading 1X40'HQ	3253CTNs(32,530Boxes)		

