



Raxwell



**RX 301 Blue Nitrile
Chemo Exam Gloves 4 mil**



Raxwell



Features

Powder Free
 No Latex
 Textured
 100 Gloves/Box
 Heavy-duty
 Fentanyl Protection
 Chemical & Oil Resistance
 Multi Purpose
 Top Quality

Standards

(US) FDA 510K / ASTM D6319/ ASTM D6978-05(2019)/ ASTM F2671-15

(EU) CE/ EN 455 part 1,2,3 &5/ EN 374 part 1, 2, 3&4/ EN 16523-1:2015/ EN ISO 21420-2020

(CN) GB4806.11-2016

Particular	Specifications	ASTM
Palm width (mm)	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X-Large	120 ± 5
Length (mm)	230 ± 5	min. 220-230
Thickness (Single Wall)	4 mil	2 mil
Tensile Strength	Before aging (MPa)	min. 21
	After aging (MPa)	min. 19
Elongation	Before aging (%)	500
	After aging (%)	500

Part No.	Size	G.W. /Carton (LBS)	Cartons/Pallet
RX3011	S	8.16	115
RX3012	M	9.04	115
RX3013	L	10.36	115
RX3014	XL	11.46	115

Packing	100Pcs/COLOR Box By Weight, 10Boxes/CTN
Box Dimension (Inch)	9.06*2.56*4.72in
Carton Dimension (Inch)	13.19*9.84*9.45in
Contain Loading 1X40'HQ	3383CTNs(33,830Boxes)



Food Service



Sanitation



Janitorial



Tattoo



Veterinary



Plumbing



Health Care



Painting



Automotive



Chemical



Industrial

EN 455




POWDER FREE



AMBIDEXTRIOUS



TEXTURED



EX Exam Grade

www.gdimed.com