



Raxwell



**RX PRO Black
Nitrile Gloves 4.5 mil**



Raxwell



Features

- Powder Free
- No Latex
- Textured
- 100 Gloves/Box
- Heavy-duty
- Fentanyl Protection
- Chemical & Oil Resistance
- Multi Purpose
- Top Quality

Standards

(US) FDA 510K / ASTM D6319/ ASTM D6978-05(2019)/ ASTM F2671-15

(EU) CE/ EN 455 part 1,2,3 &5/ EN 374 part 1, 2, 3&4/ EN 16523-1:2015/ EN ISO 21420-2020

(CN) GB4806.11-2016

Particular	Specifications	ASTM
Palm width (mm)	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X-Large	120 ± 5
Length (mm)	230 ± 5	min. 220-230
Thickness (Single Wall)	4.5 mil	2 mil
Tensile Strength	Before aging (MPa)	min. 21
	After aging (MPa)	min. 19
Elongation	Before aging (%)	500
	After aging (%)	500

Part No.	Size	G.W. /Carton (LBS)	Cartons/Pallet
RX120223-B	S	8.60	120
RX120224-B	M	9.48	120
RX120225-B	L	10.80	120
RX120226-B	XL	11.90	120
RX120241-B	XXL	13.01	120

Packing	100Pcs/COLOR Box By Weight, 10Boxes/CTN
Box Dimension (Inch)	9.06*2.36*4.53\23*6*11cm
Carton Dimension (Inch)	12.20*9.45*9.45\31x23x24cm
Contain Loading 1X40'HQ	3976CTNs(39,760Boxes)



Food Service



Sanitation



Janitorial



Tattoo



Veterinary



Plumbing



Health Care



Painting



Automotive



Chemical



Industrial

EN 455




POWDER FREE



AMBIDEXTRIOUS



TEXTURED



EX Exam Grade

www.gdimed.com