



# Raxwell



---

**RX PRO Black  
Nitrile Gloves 5.5 mil**

---



# Raxwell



## Features

Powder-Free  
 Latex-Free  
 Fully-Textured  
 100 Gloves/Box  
 Heavy-duty  
 Fentanyl Protection  
 Chemical & Oil Resistance  
 Multi-Purpose  
 Top Quality

### Standards


(US) FDA 510K / ASTM D6319 / ASTM D6978-05(2019) / ASTM F2671-15

(EU) CE / EN 455 part 1, 2, 3&5 / EN 374 part 1, 2, 3&4 / EN 16523-1:2015 / EN ISO 21420-2020

(CN) GB4806.11-2016

Particular	Specifications	
Palm width (mm)	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X-Large	≥ 115
Length (mm)		≥ 230
Thickness (Single Wall)		5.5 mil
Tensile Strength	Before aging (MPa)	min. 14
	After aging (MPa)	min. 14
Elongation	Before aging (%)	500
	After aging (%)	400

Part No.	Size	Carton Gross Weight (lbs)	Cartons/Pallet
RXF20245-B	S	10.98	105
RXF20246-B	M	11.93	105
RXF20247-B	L	13.07	105
RXF20248-B	XL	14.11	105
RXF20249-B	XXL	13.14	105
Packing		S/M/L/XL:100Pcs/Color Box By Weight, 10Boxes/CTN	
		XXL:90Pcs/Color Box By Weight, 10Boxes/CTN	
Box Dimensions			9.45x4.53x2.56in
Carton Dimensions			13.19x9.84x9.45in
Container Capacity 40' HQ		3200CTNs(32,000Boxes)	



# www.gdimed.com