



Raxwell



**RX PRO Black
Nitrile Gloves 6.5 mil**



- Powder Free
- No Latex
- Textured
- 100 Gloves/Box
- Heavy-duty
- Fentanyl Protection
- Chemical & Oil Resistance
- Multi Purpose
- Top Quality

Standards

(US) FDA 510K / ASTM D6319/ ASTM D6978-05(2019)/ ASTM F2671-15

(EU) CE/ EN 455 part 1,2,3 &5/ EN 374 part 1, 2, 3&4/ EN 16523-1:2015/ EN ISO 21420-2020

(CN) GB4806.11-2016

Particular		Specifications
Palm width (mm)	X-Small	75 ± 5
	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X-Large	115 ≥
Length (mm)		230 ≥
Thickness (Single Wall)		6.5 mil
Tensile Strength	Before aging (MPa)	min. 14
	After aging (MPa)	min. 14
Elongation	Before aging (%)	500
	After aging (%)	400

Part No.	Size	(LBS) /Carton	Cartons/Pallet
RXE20236-B	XS	11.46LBS	84
RXE20231-B	S	12.79LBS	84
RXE20232-B	M	14.11LBS	84
RXE20233-B	L	15.65LBS	84
RXE20234-B	XL	16.98LBS	84
RXE20235-B	XXL	15.21LBS	84
Packing	100Pcs/Color Box By Weight, 10Boxes/CTN		
Box Dimension	9.45X2.76x4.92in		
Carton Dimension	14.17X10.24X9.84in		
Contain Loading 1X40'HQ	3000CTNs(30,000Boxes)		

