

Raxwell



**8.5mil 12inch
Extended Cuff**

Raxwell

Features



Powder Free
 No Latex
 Textured
 50 Gloves/Box
 Heavy-duty
 Fentanyl Protection
 Chemical & Oil Resistance
 Multi Purpose
 Top Quality

Standards

(US) FDA 510K / ASTM D6319/ ASTM D6978-05(2019)/ ASTM F2671-15

(EU) CE/ EN 455 part 1,2,3 &5/ EN 374 part 1, 2, 3&4/ EN 16523-1:2015/ EN ISO 21420-2020

(CN) GB4806.11-2016

Particular	Specifications	ASTM
Palm width (mm)	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X-Large	120 ± 5
Length (mm)	300 ± 5	min. 300-305
Thickness (Single Wall)	8.5 mil	8.0 mil
Tensile Strength	Before aging (MPa)	min. 21
	After aging (MPa)	min. 19
Elongation	Before aging (%)	500
	After aging (%)	500

Part No.	Size	LBS /Carton	Cartons/Pallet
RX40001	S	4.2kgs	84
RX40002	M	4.7kgs	84
RX40003	L	5.2kgs	84
RX40004	XL	5.7kgs	84
RX40005	XXL	6.2kgs	84
Packing	50Pcs/COLOR Box By Weight, 10Boxes/CTN		
Box Dimension	11.81X1.97X4.92in/30X5X12.5CM		
Carton Dimension	10.24X10.24X12.2in/26X26X31CM		
Contain Loading 1X40'HQ	3253CTNs(32,530Boxes)		



Food Service



Sanitation



Janitorial



Tattoo



Veterinary



Plumbing



Health Care



Painting



Automotive



Chemical



Industrial

EN 455




POWDER FREE



AMBIDEXTRIOUS



TEXTURED



-25°C min

+70°C max



EX
Exam Grade

www.gdimed.com