



# Raxwell



**RX PRO Orange  
Nitrile Gloves 8.5MIL**

# Raxwell

## Features



Powder Free  
 No Latex  
 Textured  
 100 Gloves/Box  
 Heavy-duty  
 Fentanyl Protection  
 Chemical & Oil Resistance  
 Multi Purpose  
 Top Quality

### Standards

(US) FDA 510K / ASTM D6319/ ASTM D6978-05(2019)/ ASTM F2671-15

(EU) CE/ EN 455 part 1,2,3 &5/ EN 374 part 1, 2, 3&4/ EN 16523-1:2015/ EN ISO 21420-2020


(CN) GB4806.11-2016

Particular	Specifications	ASTM
Palm width (mm)	Small	85 ± 5
	Medium	95 ± 5
	Large	105 ± 5
	X-Large	115 ± 5
	2X-Large	120 ± 5
Length (mm)	230 ± 5	min. 220-230
Thickness (Single Wall)	8.5mil	8 mil
Tensile Strength	Before aging (MPa)	min. 21
	After aging (MPa)	min. 19
Elongation	Before aging (%)	500
	After aging (%)	500

Part No.	Size	LBS /Carton	Cartons/Pallet
RX50001	S	18.08LBS	56
RX50002	M	21.38LBS	56
RX50003	L	22.27LBS	56
RX50004	XL	23.37LBS	56
RX50005	XXL	24.52LBS	56
Packing	100Pcs/Color Box By Weight, 10Boxes/CTN		
Box Dimension	9.25X3.54X4.90in		
Carton Dimension	18.11X10.23X9.64in		
Contain Loading 1X40'HQ	2200CTNs(22000Boxes)		



EN 455



EX  
Exam Grade

[www.gdimd.com](http://www.gdimd.com)