



Innovative Life Sciences Tools

 major
science

2019-2020 BROCHURE

Power Supply



Power Supply

NEW



MP-310 / MP-320 / MP-510



Cat. No.	MP-100	MINI-300	MINI-500	MP-310	MP-320	MP-510
			Upgraded			
Display	N/A	LED		2.4" TFT LCD		
Output voltage / Inc.	50V/100V	10-300V/1V	10-500V/1V	5-300V/1V		5-500V/1V
Output current / Inc.	400mA	10-400mA/1mA	1-400mA/1mA	1-700mA/1mA	10-3000mA/10mA	1-800mA/10mA
Output power / Inc.	40W	60W	120W	Output: 150W Input: 200W	300W	
Output type	Constant voltage	Constant voltage/Current		Constant voltage/Current /Power		
Timer	N/A	999 (min) with alarm/Continuous		Constant mode: 9999 (min) with alarm/Continuous Programmable mode: 999 (min) with alarm/Continuous		
Rated voltage	100V-120V~ Only			100V-240V~		
Program	N/A	2-steps (V, mA, ⌚)		Pre setting: Up to 6-steps (V, mA, W, ⌚), 30 programmed files		

Nucleic Acid / Protein Electrophoresis and Blotting

NEW

Wide / Mini Horizontal Gel Electrophoresis System

NEW

CE



Cat. No	MJ-105A	MT-108
Dimension(mm)(WxLxH)	148.3x145.8x59.1	136x188x58
Gel dimension(mm)(WxL)	107x60	105x83
	52x60	50x83
Maximum sample	25 samples	
Rapid casting gel	Use gel maker stand	

* For more information, please refer to p9



www.majorsci.com

info@majorsci.com



US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107



Taiwan Office

Head Office:

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 33061, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Factory:

No.19, Ln. 207, Huakang St. Bade Dist.,
Taoyuan City 33464 Taiwan
T/ +886-3-3623319
F/ +886-3-3623133



Shanghai Office

Room 612, International business exhibition center,
9300 Hunan Road, Pudong, Shanghai, China
National toll-free No.:400-823-9177
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D. No. 12-13-99, Satguru Apts. Extn.
Street. No. 3, Lane No.1, Tarnaka
Secunderabad – 500 017. India
T/ +91-40-27001515
T/ +91-40-27001586