



PRODUCT	EMPTY FPLC 8ml	EMPTY FPLC 30ml	EMPTY FPLC 45ml	EMPTY FPLC 80ml
CAT. No.	FPLC8-3	FPLC30-2	FPLC45-2	FPLC80-1
CONTAINS	3 columns, 9 frits & 6 caps	2 columns, 6 frits & 4 caps	2 columns, 6 frits & 4 caps	1 column, 3 frits & 2 caps
CAPACITY	8	30	45	80
(ml of packed resin) <sup>(1)</sup>				
DIMENSIONS	~12 x 70	~21 x 87	~21x 137	~26 x 144
Diameter x height (mm)				
COLUMN MATERIAL	Polypropylene			
FRIT MATERIAL	Polyethylene			
FRIT PORE SIZE	10μm			
MAX PRESSURE	200/14/1.38			
_(psi/bar/MPa)				
CHEMICAL STABILITY	Stable in all commonly used reagents			
CONNECTIONS <sup>(2)</sup>	Standard connection usually supplied with the common FPLC instruments (such as ÄKTA™)			
ACCESSORIES	FPLC LUER lock-male/ CAT.No. FPLCM-1			

FPLC LUER lock-female/CAT. No. FPLCF-1

ABT TS EmptyFPLC 2017/A

<sup>(1)</sup>Packed volume depends on the type of resin — only use as a guide.
(2)For fitting the column to the FPLC equipment. If the user does not have standard connections. ABT supplies them CAT. No. FPLCM-1 /FPLCF-1 Larger amounts available on request.