

| | | <u>DNA Stain Type</u> | |
|--------------------------|---------------------------------|---------------------------------|---------------------------------|
| | | Ethidium Bromide | SmartGlow™ SYBR® Green |
| <u>Excitation Source</u> | SmartDoc blue illumination base | Not recommended | E5001-ORANGE Orange filter |
| | Blue light transilluminator | E5001-ORANGE Orange PMMA Filter | E5001-ORANGE Orange PMMA Filter |
| | UV transilluminator | E5001-590 590nm Filter | E5001-535 535nm Filter |

For optimum image quality:

- Remove the phone's case
- Line up the camera lens with the photo filter

Place the photo filter into the recess on the top platform

Remove the extender for close ups or small gels



E5001-ORANGE Orange PMMA Filter (use with blue light)



E5001-535 535nm Filter (use with UV)



E5001-590 535nm Filter (use with UV)