

Caremed Alrick.

Wireless fall prevention

Alerta manufacture wireless fall prevention sensors, delivering solutions to proactively maximise user safety whilst not compromising on comfort



WIRELESS DELUXE ALERTAMAT



A robust luxury carpet mat, most effective when placed at the bedside or doorway; activating the alarm when a user leaves their safe zone by stepping on the mat.

- ✓ Selecting one of 5 colours so that it blends with the existing flooring makes it an ideal solution for situations where a dementia user may otherwise try to avoid the mat.
- ✓ Can be connected to a nurse call system for centralised activation, or used standalone.
- ✓ X5 colours to choose from.
- ✓ Includes mat and wireless transmitter.

*Wall point receiver sold separately.

WIRELESS FLOOR+ ALERTAMAT



This beige coloured Alerta floormat is most effective when placed by the bedside or at a room doorway, thereby activating the alarm when a user leaves by stepping on the mat. The large size means it is hard to avoid by over-stepping.

- ✓ Wireless systems can be connected to a nurse call system for centralised activation, or used on their own, independent of a nurse call system to give localised alerts.
- ✓ Dementia friendly.
- ✓ Dimensions: 60 x 90 cm (23.5 x 35.5").
- ✓ Includes mat and wireless transmitter.

*Wall point receiver sold separately.

WIRELESS BED ALERTAMAT



The Bed Alertamat is most effective when concealed under the bed mattress; activating the alarm when a user leaves their bed, and resetting automatically when the user returns to the bed.

- ✓ Wireless systems can be connected to a nurse call system for centralised activation, or used on their own, independent of a nurse call system to give localised alerts.
- ✓ Hygienic, antibacterial easy-clean material.
- ✓ Dimensions: 25 x 76 cm (10 x 30").
- ✓ No cables, extra safe and versatile.
- ✓ Includes mat and wireless transmitter.

*Wall point receiver sold separately.

WIRELESS CHAIR ALERTAMAT



The Chair Alertamat is most effective when hidden under the chair cushion; activating the alarm when a user leaves their chair or slowly sinks down in the chair. An alert is given to the carer to check on the resident. The alarm is reset automatically when sits back into the chair.

- ✓ Wireless systems can be connected to a nurse call system for centralised activation, or used on their own, independent of a nurse call system to give localised alerts.
- ✓ Hygienic, antibacterial easy-clean material.
- ✓ Dimensions: 25 x 38 cm (10 x 15").
- ✓ Wireless connect-and-play system.
- ✓ Includes mat and wireless transmitter.

*Wall point receiver sold separately.

WIRELESS CRASH MAT, INC WIRELESS RECEIVER WIRELESS WALL POINT RECEIVER



Designed to be placed next to the bed to minimise injury if the user should fall from bed and alert the carer, the mat incorporates sensor pads and a wireless transmitter within the foam which connects to the Alerta wall point receiver.

- ✓ Includes wireless receiver.
- ✓ It can be connected to a nurse call system for centralised activation, or used on its own, independent of a nurse call system.
- ✓ The volume can be deactivated for silent in-room operation.
- ✓ Helps prevent a user remaining on the mat undetected for an unknown length of time.
- ✓ PU wipe clean surface.
- ✓ Unfolded the mat is 180cm long x 60cm wide x 6cm high.

WIRELESS ALERTA NURSE CALL BUTTON



The call button provides the resident with reassurance that help is at hand. With a large easy to press button it activates the alarm when a user presses the button. This alerts the carer to check on the resident.

- ✓ Lanyard and wall mounted holder included.
- ✓ Wireless connect-and-play system.
- ✓ Takes 1x12v 23A battery.

*Wall point receiver sold separately.



Can be connected to a nurse call system for centralised activation or used on its own, independent of a nurse call system to give localised alerts through this device only. The volume can be changed for silent in-room operation when connected to a nurse call system.

- ✓ Can be wirelessly paired with up to 5 Alerta wireless products, please view the user manual [here](#).
- ✓ It uses radio frequency technology to connect.
- ✓ Wireless range of up to 95 metres.
- ✓ Wireless connect-and-play system.
- ✓ High, low and silent volume control.
- ✓ Low battery alert.
- ✓ Takes 2x1.5v AA batteries.

WIRELESS ALERTA DETECT



The Alerta Detect uses infrared technology to activate the alarm when body movement is detected within the range of the infrared beam. The Alerta Detect is most effective when mounted at the door, inside a bathroom (if this is a typical place for falls to take place) or by a bedside; activating the alarm when a resident leaves their position.

- ✓ Can be connected to a nurse call system through the Alerta Wall Point Receiver for centralised activation, or used on its own, independent of a nurse call system.
- ✓ Wireless connect-and-play system.
- ✓ No cables, extra versatile.

*Wall point receiver sold separately.