

## Manufacturer's Specifications

Brand	KANBO
Model Name	KN95 Folding face mask (Non-Medical)
Certification	CE-EN149:2001+A1:2009 / FDA#10063729 /GB2626-2006
Layers	4 layers
Color	White
Function	Anti-Pollination/ anti-dust/ anti-disease/ anti-saliva
Construction	Surface <u>spunbond</u> + middle <u>meltblown</u> + hot air cotton+white <u>spunbonded</u> cloth + ear band + nose bridge
Mask Size	W108 X H160mm ( <u>Earloop</u> )
Packing	1pcs/bag, 50bags/box, 20boxes/carton
Box size/ Weight	L27*W19*H10.7cm; Weight: 440G
Carton size/ Weight	L56.5*W40*H60cm; Weight: 9.8KG

# User Manual



**KANBO®**

**KN95**

**FOLDING FACE MASK**

SELF-PRIMING FILTERANTH-PARTICLE RESPIRATOR

Bacterial filtration efficiency (bfe) ≥95%

**PRODUCT DESCRIPTION**

Product Name	KN95 Folding face mask (Earloop)
Packing Size	L270mmxW190mmxH107mm
Color	White
Main Function	Respiratory protection for non-oil dust environment
Storage Period	2 years

**STORAGE CONDITION:** Unwashed products under conditions of guaranteed storage temperature of 20°C to 30°C and relative humidity of less than 80%.

**DATE OF PRODUCTION:** On package

**BATCH NUMBER:** On package

**IMPLEMENTATION CRITERIA:** GB2626-2006 EN149:2001

CE ISO system certified company **1** year

**⚠ WARNINGS**

- [1]** It is recommended that masks be replaced every 3 days when the air quality index is Grade 3 (high pollution) and below. It is recommended to replace the mask every two days when the air quality index is level 4 (moderate pollution) and below. It is recommended to replace the mask every day when the air quality is level 5 (heavy pollution) and below. To achieve better protection, please follow the correct usage instructions.
- [2]** Can not be used in toxic gas environment (can not replace gas masks), only for non-oil particulate matter respiratory protection.
- [3]** Do not use this mask when the pollutant concentration reaches immediate life-threatening or health-related concentration; leave the contaminated area immediately in case of breathing difficulties or dizziness or other discomfort.
- [4]** New masks should be replaced when they are damaged or when breathing resistance becomes excessive.
- [5]** Shall not modify, clean, abuse or misuse this mask.
- [6]** Found moldy, damaged packaging, after expired strictly prohibited to use.

**PRODUCTION ENTERPRISE**  
Shenzhen Horb Technology Corp. Ltd.

**ADDRESS**  
Block B, Faming 16-tech Industrial Park, Xintian Rd,  
Fuyong Town, Bao'an District, Shenzhen, China.

**POST CODE:** 518103

**WEBSITE:** www.horb.com.cn

**TEL:** 0086-755-29491987 29236367

**EMAIL:** sales@horb.com.cn

MADE IN CHINA

## KN95 Manual

Size	W65mmxH100mm
Material	Book paper
Printing color	Single color printing
Product code	5.10.04.0001

# FDA Certificate



## CERTIFICATION OF REGISTRATION

**This certifies that:**

**Shenzhen Horb Technology Corp., Ltd.**

**4-7 Floor, B Block, Funing Hi-Tech Industrial Park, No.71-2 Xintian Rd, Fuyong  
Town, Bao'an District, Shenzhen, Guangdong, CN 518000**

is registered and has listed the following medical device with the U.S. Food and Drug Administration :

**Owner/Operator Number : 10063729**

**Listing Number:** [REDACTED]

**Product Code : MSH**

**Product : Respirator, surgical - Folding Protective Mask**

**Model(s) : N95, KN95**

**Date Of Registration Status: 2020**

Registrar Corp will confirm that such registration remains effective upon request and presentation of this certificate until the end of the year stated above, unless said registration is terminated after issuance of this certificate. Registrar Corp makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. Registrar Corp assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding."

The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. Registrar Corp is not affiliated with the U.S. Food and Drug Administration

  
Chief Engineer / Zhaohao Lee

Issued: Mar. 23, 2020

Validity Period: 2020-12-31



Web rep/www.fda.gov Tel 1-888-INFO-FDA(1-888-0332) E-mal vcrp@fda.gov

# CE Certificate



# Certificate

## ICR Polska/CE/V/00RE256

**Name and address of certificate owner:** Shenzhen Horb Technology Corp., Ltd.  
Floor 4-7, Block B, Funing Hi-Tech Industrial Park, No.71-2 Xintian Rd,  
Fuyong Town, Bao'an District, Shenzhen City, 518000, China.

**Name and address of manufacturer:** Shenzhen Horb Technology Corp. Ltd.  
Floor 5-7, Block B, Funing Hi-Tech Industrial Park, No.71-2 Xintian Rd,  
Fuyong Town, Bao'an District, Shenzhen City, 518000, China.

**Product name:** KN95 Mask

**Product types:** 1.7.02.02.0001

**Product trademark:** n/a

- This certificate confirms that the product meets the requirements of the following standard
- The conformity was demonstrated based on following standard(s) listed by European Commission as harmonized with **Regulation/Directive 2016/425**

EN 149:2001+A1:2009

The certification has been carried out in accordance with individual rules and conditions agreed with the applicant. Evaluation has been carried out in accordance with:

**Test report(s) No.:** CTO200312018LRS

**Test conducted by:** Shenzhen CTO Technology Service Co., Ltd

**Certificate issue date:** 02.04.2020

**Expiration date:** 01.04.2025

NOTE:

- This certificate refers to the above mentioned product and its conformity in regards of above mentioned standard(s) was proven on test sample
- This certificate does not imply meeting all essential requirements, assessment of the series-production or any other restricted Notified Bodies conformity assessment procedure appropriate for the product
- This certificate holder shall use this certificate in connection to declaration of conformity and technical data relevant for the product the certificate was issued

Certificate modifications: n/a



**ICR Polska Co. Ltd.**  
Plac Przymierza 6  
03-944 Warsaw  
www.icrpolska.com  
e-mail: [icrpolska@icrqa.com](mailto:icrpolska@icrqa.com)

Director: Rafał Kalinowski

Warsaw, 02.04.2020



# GB2626-2006 Certificate



## 检验检测报告 (电子版)

报告查询网址: [www.gtgc.net.cn](http://www.gtgc.net.cn)

防伪码: MJX-6718-24

立即下载



No.: 200073437

共3页 第1页

委托单位	深圳市中明科技股份有限公司 地址: 深圳市宝安区福永街道新田大道71-2号立群科技园高新产业园5栋4-1层		
客户认定信息	D/94 防雾式口罩 30个 商标: KAISO 型号: S.T.02.02.0001		
检验性质	委托检测	样品受理/测试开始日期	2020-04-14
		报告签发日期	2020-04-20
判定依据	GB 2626-2006 《呼吸防护用品 自吸过滤式防颗粒物呼吸器》		
综合检验结论	---		
检验检测结果	检验检测项目	判定依据	判定
	NaCl 颗粒物过滤效率	GB 2626-2006	符合
	吸气阻力	GB 2626-2006	符合
	呼气阻力	GB 2626-2006	符合
	可燃性	GB 2626-2006	符合
	头带	GB 2626-2006	符合
备注	本报告中检验检测项目均在相应标准规定的环境条件下进行(有注明的除外)。复印件、副本未经加盖报告书确认章无效。本报告检验检测地址为: 广州市番禺区珠江路1号。		



签发: 郭子山 高级工程师 *郭子山*

总部: 广州市番禺区珠江路1号  
花都实验室: 广州市花都区梯柯镇梯柯村西岗1号

电话: 020-61994598/01994599  
电话: 020-37721161

## Qualified Certificate

**KANBO®**

**合格证**

**QUALIFIED CERIFICATE**

生产企业: 深圳市中明科技股份有限公司  
地址: 深圳市宝安区福永街道新田大道福宁高新产业园B栋  
Production enterprise:  
Shenzhen Horb Technology Corp.Ltd.  
Address: Block B, Funing Hi-Tech Industrial Park, Xintian  
Rd, Fuyong Town, Bao'an District, Shenzhen, China.

**KN95折叠式口罩(非医用)**  
**Kn95 folding face mask(non-medical use)**

品牌/Brand: KANBO

型号/Model: 1.7.02.02.0001

规格/ Specification: 108mm\*160mm

数量/Quantity: 50 pcs

检验员/Inspector:

有效期/Valid Period: 2 Years

执行标准/Standard:

GB2626-2006 EN149-2001+A1:2009

主要成分/Main ingredients:

无纺布、熔喷布、热风棉/non-woven fabric,  
meltblown cloth, hot air cotton



-----  
生产日期/Manufacturing Date:

生产批号/Batch No.: