


Classic Spinal Orthosis Spina™ IV Model 1400

Easy to fit, high patient acceptance
Approved Code L0637

Traditional style spinal orthosis lined with a soft-foam removable liner. Apron front is soft-foam material with pocket for moldable insert. Orthosis attaches with 2" double sided Velcro webbing – straps are adjustable at both ends.

Standard – 14" back height & 9" front
Options – low profile front & pendulous insert

Size	Hip Circumference
Small.....	30" – 35"
Medium.....	36" – 41"
Large.....	42" – 47"
Ex-Large.....	48" – 52"



 **U.S. Orthotics**®
INC.

800-825-5228

8605 Palm River Rd.
Tampa, FL 33619
Ph: 813-621-7797
Fax: 813-623-1055
Toll-Free: 800-825-5228

www.usorthotics.com

Night Splints
Elastic Back Supports
Neoprene Supports
Knee Immobilizers
Walking Boots

orders@usorthotics.com

AFO
Maternity Belts
Ankle Supports
Wrist Supports
Side-Lace Corsets

info@usorthotics.com

Pediatric Supports
Clavicle Straps
Arm Slings
Cervical Collars
Rib Belts & Binders
Prosthetic Suspension

Spina™ Kydex Chairback
Spina™ II Plastic Orthoses
Hyperextension Orthoses
Pre-Form Foot Orthotics
Shoulder Immobilizers