

V-Line[®] Spinal Orthoses



Our professional grade V-Line[®] Spinal Orthoses series features comfort, stability and a wide array of options. A single piece design makes for easy application and effective compression. Pre-molded posterior inserts come in 0° – 7° – 12° – 25° lordosis with a spinal relief cut-out. Pre-molded anterior inserts comes neutral or pendulus and are pre-slotted to accommodate an anterior thoracic extension. The anterior pocket is fully adjustable providing correct positioning and accommodating a wider range of sizes. Designed with two smooth cross-pull compression straps, with elastic tension ends. Hand-pull tabs on all ends.



PRO-V1 **L0637**
14" back height – LSO
Standard 9" front – Model V1
Low profile 7" front – Model V2



PRO-V3 **L0627**
12" back height – LSO
Standard 7" front




PRO-V9 **L0631**
14" back height – LSO
Standard 9" front – Model V9
Low profile 7" front – Model V10
2 piece design



PRO-V6 **L0456**
18" back height – TLSO
Standard 9" front – Model V6
Low profile 7" front – Model V7

Anterior panel options
Neutral or Pendulous



Posterior panel options
0° – 7° – 12° – 25° lordosis



Dorsolumbar Extension-DLE
Fits all V-Line[®] PRO-V style LSO spinal orthoses



Anterior Thoracic Extension-ATE
6" and 8" sternal shields



PRO-V20 **L0456**
Model V6 with the addition of an ATE

SIZE	HIPS
Small	32" – 38"
Medium	36" – 42"
Large	40" – 46"
Ex-Lg	44" – 50"
2XL	48" – 54"
3XL	52" – 58"
4XL	56" – 62"