

## Section 1: Product and company information

### Product identifier

<b>Product name</b>	<b>Chitozen 5-coverslip Standard Kit</b>
<b>Substance name</b>	N/A
<b>Reach No.</b>	A registration number is not available for this substance as the substance or its uses are exempted from registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline

<b>HTS Code</b>	3822.00
-----------------	---------

### Relevant identified uses of the substances or mixture and uses advised against

#### Relevant Identified uses

<b>Main use category</b>	Laboratory use
<b>Professional use spec</b>	R&D use only
<b>Identified used</b>	Laboratory chemicals

### Details of the supplier of the safety data sheet

<b>Company</b>	Idylle 68 boulevard de Port-Royal 75005 Paris France
<b>Technical phone</b>	+33184172227
<b>Email address</b>	<a href="mailto:contact@idylle-labs.com">contact@idylle-labs.com</a>


### Emergency phone number

<b>France</b>	112
<b>Other countries</b>	ORFILA (INRS) : + 33 (0)1 45 42 59 59

## Section 2: Technical information

### Kit and product description

- Each kit is composed of **5 chitosan-coated coverslips** and packaged in a shockproof box in which up to 5 standard units can be stored.
- Each Chitozen Standard coverslip is composed of the **chitosan-coated coverslip**, standing against a **protective slide** on which the product label stands. The chitosan-coated coverslip is originally linked to the protective slide by two removable adhesives. These two removable adhesives are to be peeled off prior to use (see protocol). The chitosan-coated side of the coverslip is the one facing up (see scheme below).
- Dimensions of the chitosan-coated coverslip (w x l x h): 75 x 25 x 0.175 mm<sup>3</sup>



### Shipping, storage and expiration

- Each kit of **5 chitosan-coated coverslips** is shipped at room temperature and shielded from direct sunlight light in the shockproof box.
- The **chitosan-coated coverslips** should be stored at room temperature, away from moisture and shielded from direct sunlight, in the shockproof box they are originally packed.
- The **chitosan-coated coverslips** can be stored and used up to 12 months after production. Production year and month is indicated by the batch code

*According to Regulation (EC) No. 1907/2006(REACH) with its amendment Regulation (EC) No. 453/2010-830/2015*

YYMM (where YY stands for the last two digits of the year and MM for the month).