



Tangible Hydra-PEG Compatible Solutions List



OUS Compatible Solutions

Cleaning solutions that are distributed outside of the United States, and compatible with Tangible Hydra-PEG.

Manufacturer	Solution Name	Solution Type
Abatron	Quattro	Multipurpose
Abatron	Amclair Weekly Enzymatic Cleaning Tablets	Enzymatic Cleaning Tablets
Abbott Medical Optics/ Johnson & Johnson	Total Care Cleaner	Cleaner
Abbott Medical Optics/ Johnson & Johnson	Total Care Wetting	Storage and Wetting Solution
Agepha Pharma	Viscosan	Eye Drops and Rewetting Drops for RGP
Alcon	AO Sept Plus	Peroxide
ASDA	All in One	Multipurpose
Avizor	Avizor All Clean Soft	Multipurpose
Avizor	Avizor Ever Clean	Peroxide and Neutralization tablets
Avizor	Avizor Ever Clean Plus	Peroxide and Neutralization tablets
Avizor	GP Multi	Multipurpose
Bausch + Lomb	EasySept Hydro+	Peroxide
Bilosa	BiloPlus	Multipurpose
Boots Pharmaceuticals	Boots Total Solution	Multipurpose
Contopharma	Contopharma Tabinone	Peroxide and Neutralization Tablets
Contopharma	Contopharma Conventional Cleaner	Cleaner
Contopharma	Contopharma GPHCL Cleaning C: Cleaning Solution for RGP Contopharma GPHCL Rinsing R: Rinsing and storage solution for RGP	Contact Lens Care System: Cleaner, and Rinsing/Storage Solution

Manufacturer	Solution Name	Solution Type
Contopharma	Conventional Rinsing	Rinsing Solution
Consol	Extend	Multipurpose
Coopervision	Coopervision Refine One Step	Peroxide
Consol	Extend All-In-One Contact Lens Solution	Multipurpose
Disop	Hidro Health	Peroxide and Neutralization Tablets
Hecht	Concare Enzyme tablets	Enzymatic Cleaning Tablets
Hecht	Sichtbar Besser	Storage and Wetting Solution
HyLo	HyLo Care	
HyLo	HyLo Comod	
HyLo	Hylo- Fresh	Multipurpose
Lobob	Reiniger	Cleaner
Lobob	Conservation	Storage and Wetting Solution
Menicon	Menicare Plus	Multipurpose
Menicon	Menicare Pure	Peroxide
MPG&E	Perfect Aqua Plus	Peroxide
Ophtecs	Cleadew MPS Advanced Care System for Soft	Multipurpose
Ophtecs	Cleadew Advanced Care System for Soft	Disinfecting Solution and Enzymatic Cleaning Tablets
Ophtecs	Cleadew Advanced Care System for GP	Disinfecting Solution and Enzymatic Cleaning Tablets
Ophtecs	Cleadew SL	Multipurpose
Oté	Oté Clean Multipurpose Solution	Multipurpose

Manufacturer	Solution Name	Solution Type
Oté	Oté Monosept+ 12	Peroxide
Prolens	Topocare GP Sensitive	Multipurpose
Protecta Vision	Ottica Napoli	Multipurpose
Rubicon	Clears	Multipurpose
Rubicon	Clear Sol	Rinsing and Fitting Solution
Sauflon	Delta Plus	Storage and Wetting Solution
Sauflon	Delta Daily Cleaner	Multipurpose
Sauflon	Proclean	Cleaner
Vita Research	Regard K	Multipurpose

*Generic Peroxide Systems may be used or recommended

Current List Of Compatible Rewetting Drops

Any drop indicated for use with contact lenses, or is preservative-free can be used and recommended.

Current List Of Chemically-Compatible Solutions with IFU Issues

Cleaning Solutions "Chemically Compatible" with Tangible Hydra-PEG but are recommended to use with incompatible products or tap water. These products will not damage the coating on their own.

Manufacturer	Solution Name	Solution Type	Notes
Bausch + Lomb	Boston One Step Liquid Enzymatic Cleaner	Liquid Enzymatic Cleaner	Instructions recommend use with Boston Advance Cleaner, which will destroy the coating

Current List Of Non-Compatible Solutions

Cleaning solutions that are not-compatible with Tangible Hydra-PEG. Do not use or recommend:

- Any alcohol based care solutions
- Any abrasive based care solutions
- Solutions listed below

Manufacturer	Solution Name	Solution Type
Avizor	GP System Cleaner	Cleaner
Avizor	GP System Conditioner	Storage and Wetting Solution
Bausch + Lomb	Boston Advance Cleaner	Advance Formula Cleaner
Bausch + Lomb	Boston Laboratory Lens Cleaner	Laboratory Use Cleaner
Bilosa	Bilosafe	Storage and Wetting Solution
Hecht	Concare Aufbewahrungslosung	Storage and Wetting Solution
Menicon	Progent	Protein Removal Cleaner
Oté	Oté Fine 40 All-in-One	Multipurpose
Oté	Oté Perm	Storage Solution