

Ascend[®] Full Face

Ascend to a higher level of comfort and convenience

Rise to a higher level of comfort and convenience with the Ascend Full Face Mask. Our patented AIR[°]gel replaceable cushion is easily detached and reattached for cleaning and cushion replacement, saving both time and money.

Touchless Spacebar

Unique design provides superior stability, maintaining a seal without contacting the skin.

Replaceable Cushion

Now you can replace the cushion without having to replace the whole mask. Available in three sizes.

S M L

Quick Release Tab

Lightweight breathable headgear removes easily with one hand.

AIR[°]gel[®]

Maximum comfort and minimal headgear tension. AIR[°]gel[®] minimizes pressure points and skin irritation.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Active Headgear Connector

Moves as the patient moves for a reliable seal.

What makes our Ascend Full Face Mask different?

Replaceable Cushion - Available in three sizes with simple and breathable quick-release AFit headgear. Our easy-to-use replaceable cushion allows you to extend the life of your mask.

Active Headgear Connector - Moves as patient moves for a reliable seal.

Touchless Spacebar - Unique design provides superior stability, maintaining a seal without contacting the skin.



AIR°gel® Cushion Technology

Our patented AIR°gel technology with an ergonomic design reduces pressure while maintaining an effective seal.

Easy-to-Replace Cushion **S M L**



Assembly: Place the top alignment pin in the notch. Push the bottom until all pins are secure.



Disassembly: Squeeze the sides of the cushion frame to release from mask.

Ascend Full Face Masks with AFit Headgear Specifications:

Product Number	Description	Cushion Size	Dimensions
50825	Mask System & headgear	S,M,L	NA
50834	Mask System with attached headgear	S,M,L	NA
50801	Mask with attached headgear	Small	3.92"D x 3.86"W x 7.00"H
50802	Mask with attached headgear	Medium	4.16"D x 4.20"W x 7.00"H
50803	Mask with attached headgear	Large	4.24"D x 4.30"W x 7.00"H
50813	Mask without headgear attached	Small	3.92"D x 3.86"W x 7.00"H
50814	Mask without headgear attached	Medium	4.16"D x 4.20"W x 7.00"H
50815	Mask without headgear attached	Large	4.24"D x 4.30"W x 7.00"H
50207	AIR°gel cushion	Small	3.92"D x 3.86"W x 7.00"H
50208	AIR°gel cushion	Medium	4.16"D x 4.20"W x 7.00"H
50209	AIR°gel cushion	Large	4.24"D x 4.30"W x 7.00"H
55058	AFit headgear	NA	NA

Also by Sleepnet®
Ascend Nasal
 featuring breakthrough AIR°gel® technology



Contact and support

Sleepnet representative in your area

www.sleepnetmasks.com

Phone: 1-800-742-3646

The following product is protected by US patents 10,974,007, 9,199,054 and other patents pending



Does not contain natural rubber latex

