

The V3 combines AIR^ogel[®] technology with acute care performance for a revolutionary home care experience. The V3 fits quickly and easily while enhancing patient comfort.

AIR^ogel[®]

Maximum comfort, minimum weight. Our patented AIR^ogel minimizes pressure points and skin irritation.

AIR^ogel[®] Forehead Cushion

Forehead Adjustment Screw

Quick-Release Slot

Lightweight, breathable headgear slides off easily with one hand.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Oxygen Port

Allows easy delivery of supplemental oxygen.

Size Indicator

Available in four sizes

S M L XL



What makes our V3 Mask different?

- Adjustment Screw -provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the bridge of the nose.
- AIR^ogel Forehead Pad - stabilizes the mask comfortably on the patient's face during a night's sleep minimizing pressure points.
- Available in four sizes with simple and breathable quick-release headgear.

Featuring Patented AIR^ogel[®] Technology

Our patented AIR^ogel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



V3 Masks with EZfit Headgear For Comfortable & Secure Treatment at Night



V3 Product Specifications:

Product Number	Style	Leak Port	Entrianment Valve	Size	Mask Weight	Mask Dimensions
50863	Vented	✓	✓	Small	109.2g	3.60"D x 4.10"W x 7.25"H
50864	Vented	✓	✓	Medium	114.3g	3.60"D x 4.30"W x 7.40"H
50865	Vented	✓	✓	Large	119.4g	3.60"D x 4.38"W x 7.60"H
50866	Vented	✓	✓	Extra Large	120.1 g	3.60"D x 4.42"W x 7.80"H
55051	EZfit Headgear	NA	NA	One size	NA	NA



Does not contain natural rubber latex

Contact us today.

General inquiries email
info@sleepnetcorp.com or call 1 (603) 758-6600

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive,
 Hampton, New Hampshire 03842 USA



sleepnetmasks.com