# Nonny™ Pediatric Nasal Mask

## Optional Extension Straps

(Included in Small Kit ONLY)

velcro close

Inno

### A. Five Point Headgear

- · Additional support point helps stabilize headgear
- Soft and comfortable
- Breathable against the skin
- Easy to wear, adjust and clean

#### **B. Forehead Pad**

- · Easy fitting, one size pad
- · One-piece design stabilizes mask on the face

#### C. Silicone Facial Cushion

• Uniquely designed for stable, superior comfort • Available cushion sizes:

Small Kit w/Small & Med. (Includes 2-Optional Extension Straps)

#### Large Kit w/Large & X-Small

(Although it is labeled "X-Small", this is the largest of the ped. cushion sizes. Think of it as "Adult X-Small")

- Easy to attach and remove from frame
- · Comfortable & hypoallergenic

#### D. Buckle w/quick snap design

- · Quick and easy to fasten/unfasten mask
- Optional buckle can be used with ball joint or Velcro<sup>®</sup>

#### E. AG Flex™ Tube

- Standard 22 mm connection
- Flexibility improves compliance



Mask Kit Description	Kit Part #	Headgear	Replacement Headgear Part #	Replacement Cushion Part #
Nonny Large Mask Kit	AG-PEDKIT-L	Large 5-Point	AG-PEDKIT-HGL	AG-PEDKIT-CXS AG-PEDKIT-CL
Nonny Small Mask Kit	AG-PEDKIT-S	Small 5-Point	AG-PEDKIT-HGS	AG-PEDKIT-CM AG-PEDKIT-CS



M\_NonnyMask Rev. 2