



Bubble Humidifier

Item#
RES050

20/case

- 6 PSI pressure relief valve
- Flows from 2-12 LPM
- Holds 500ml
- Disposable, durable plastic
- Unique diffuser designed to provide patient with comfortable, therapeutic humidity and minimal noise
- Defined maximum and minimum water level lines
- Recessed nipple reduces possibility of breaking off outlet connection
- Holds 500ml
- Preset audible alarm at 6 PSI pressure relief
- Allows for use of longer lengths of oxygen tubing without the audible alarm sounding

Tubing Length	Recommended Max Flow Rate
Up to 20ft	12 LPM
21-30ft	10 LPM
31-40ft	9 LPM
41-50ft	8 LPM