

iGO2 Portable Oxygen Concentrator (POC)

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Conserver Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market
- 8.4" (H) x 3.5"(W) x 8.6"(L)
- Portable Oxygen Concentrator iGo family

| PRODUCTS | Description |
|------------------|--|
| 125D | iGO2 Portable oxygen concentrator NA |
| 125D-XB | iGO2 Portable oxygen conc w/xtra bat NAC |
| SPECIFICATIONS | Description |
| Battery Duration | 3.5 hours @ setting 2 |
| HCPCS | E1390 + E1392 |
| Oxygen Purity | 90% +4/-3 over entire operating range |
| Sound Level | 37db |
| Warranty | Drive DeVilbiss warrants this equipment to be free from defects in workmanship and materials for five (5) years from date of factory shipment to the original purchaser (typically the healthcare provider) with the exceptions of the Compressor (warranted for 3 years), Sieve Beds (warranted for 1 year), and Battery (warranted for 1 year). The carry case is excluded from this warranty. These warranties are limited to the Buyer of new equipment purchased directly from Drive DeVilbiss, or one of its Providers, Distributors, or Agents. Drive DeVilbiss obligation under this warranty is limited to product repair (parts and labor) at its factory or at an Authorized Service Center. Routine maintenance items, such as filters, are not covered under this warranty, nor does it cover normal wear and tear. |
| ACCESSORIES | Description |
| DV51D-606 | US POWER CORD, EA |
| DV6X-619 | DC car adaptor for IntelliPAP 2, EA |
| 125D-613 | iGO2 Battery pack, EA |
| 125D-670 | iGO2 Carrying case, EA |
| 125CH-613 | iGO2 Battery Charger (US), EA |