

SleepWeaver® Anew™ Online Sizing Template

Instructions for Use:

1. To ensure the Sizing Guide is true to scale, verify the printer is set to print at 100%.
2. After printing, check the size of the printed version of the Sizing Guide by lining up a ruler with the ruler scale on the guide.
3. Carefully cut out the Sizing Guide on the dotted line.
4. Remember, a soft fabric mask is different than a hard plastic mask. Fabric can be pulled tighter and relaxes as it is used daily.
5. If you fall "in between" sizes, please order the smaller size of Anew.

SleepWeaver®
Anew™

Sizing Guide

4 in

3 in

2 in

1 in

0 in

B : Bridge of Nose

LARGE 100965 ↗

REGULAR 100960 ↗

SMALL 100955 ↗

	Small	Regular	Large
Nasal Width	1 in 2.5 cm	1.25 in 3.2 cm	1.5 in 3.8 cm
Height	3.25 in 8.2 cm	3.5 in 8.9 cm	3.75 in 9.5 cm

A : Middle of Chin ↗

888-505-0212 DirectHomeMedical.com