

# READY HEAT

# THE READY-HEAT FIRST RESPONDER LINE

*The Temperature Management Solution*

The First Responder solution for active patient warming in any casualty situation:



Emergency Rescue




Disaster Preparedness



Air Transport


## Ready-Heat II

- The #1 selling active warming blanket
- Quick warming blanket, with 8+ hours of active warmth
- Ideal for all patient transport and emergency rescue
- Water resistant top layer made of coated medical grade non-woven material and a foam core
- Also available in Torso size




## Ready-Heat Disaster Blanket

- For use during disasters and mass casualty events
- Deploys in seconds and heats in minutes
- Supplies consistent heat for 8 hours
- 5 year shelf life
- Water resistant top layer made of coated medical grade non-woven material and a foam core



## Ready-Heat Infant Cocoon

- Designed to actively warm the most precious cargo
- Consists of 6 thermal pads strategically placed to warm infant's body
- Hood assists in the retention of body warmth and the padded foam core and dual hook and loop design provide maximum comfort and fit



## For the First Responder

Stay warm on the job with proprietary Ready-Heat technology in a disposable heated vest.

## Include Ready-Heat in Your:

- First Aid Kits
- Wilderness Survival Kits
- Mountain Rescue Kits
- Pediatric Transport Kits
- Water Rescue Kits
- Trauma Kits



MANUFACTURED BY TECHTRADE LLC

30 Montgomery Street • Suite 690 • Jersey City, NJ 07302 • 1-866-491-TECH  
Fax: 201-706-8131 • info@ready-heat.com • www.ready-heat.com



## 4 PANEL BLANKET

NSN #: 6532-01-525-4062

Part Number / Per Case: S4RHMD / Case of 8



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 34" x 48" / 86 cm x 122 cm


**Weight:** 1.75 lbs. / 794 gm

**Description:** Torso Sized Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250 lbs./115 kg). Military Medical Protocol.

## 6 PANEL BLANKET

NSN #: 6532-01-525-4063

Part Number / Per Case: S6RHLG / Case of 8



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 34" x 60" / 86 cm x 152 cm

**Weight:** 2.5 lbs. / 1.13 kg

**Description:** Full Body Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250lbs./115 kg). Number 1 choice of US Air Force transports.

## 1 PANEL BLANKET

NSN #: 6532-01-538-1525

Part Number / Per Case: S1RHSM / Case of 12



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 16" x 18" / 41 cm x 46 cm


**Weight:** 7 oz. / 198 gm

**Description:** Mini Blanket. Suitable for injured extremities. Fully adjustable with 2 hook and loop closures. Perfect for First Aid and Survival Kits.

## READY-HEAT II BLANKETS

NSN #: 6532-01-538-0759

Part Number / Per Case: G12RHZ / Case of 8\*



**Warming Temperature:**  
100°F / 38°C in 8-10 Minutes

**Maintains Heat:** 8 Hours

**Size:** 34" x 48" / 86 cm x 122 cm

**Weight:** 1.5 lbs. / 794 gm


**Description:** First Responder Blanket. Used to Treat Shock, Trauma and Pre-Hypothermia. Has a foam core for water resistance and is lightweight and comfortable. Quick to reach operating temperature. Works well at high altitudes.



\* Also available in Torso Size  
G6RH2T / Case of 12  
34" x 24"/86 cm x 61 cm, 14 oz./397 gm

## TEMPERATURE MANAGEMENT BLANKETS

Part Number / Per Case: SB9RH9120 / Case of 10\*



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 34" x 60" / 86 cm x 152 cm

**Weight:** 18 oz. / 510 gm

**Description:** Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Prevents hypothermia, increases patient comfort for higher quality of care, effective solution for infection control.



\* Also available in Half Body  
SB6RH9120 / Case of 20  
34" x 24"/86 cm x 60 cm, 12 oz./350 gm

## READY-HEAT PERSONAL USE

Part Number / Per Case: GRHV06A / Case of 12



**Warming Temperature:**  
100°F / 38°C in 10-15 Minutes

**Maintains Heat:** 8 Hours

**Size:** Overhead Style with Adjustable Hook and Loop Straps. One Size Fits All.

**Weight:** 1.75 lbs. / 794 gm

**Description:** Designed to keep the First Responder warm during long missions. Also great for outdoor activities (Hiking, Skiing, Camping, Hunting, Fishing, Marathons and Mountain Climbing).

## INFANT WARMING MATTRESS

Part Number / Per Case: S2RHIM / Case of 12



**Warming Temperature:**  
100°F / 38°C in 15-20 Minutes

**Maintains Heat:** 8 Hours

**Size:** 11" x 13" / 28 cm x 33 cm

**Weight:** 7 oz. / 198 gm

**Description:** Infant Temperature Management. Used in conjunction with the Infant Warming Cocoon for constant warming. Great for infant transport within hospitals.

## INFANT WARMING COCOON

Part Number / Per Case: GIW6C / Case of 12



**Warming Temperature:**  
100°F / 38°C in 15 Minutes

**Maintains Heat:** 8 Hours


**Size:** 31.5" x 28" / 80 cm x 71 cm

**Weight:** 12 oz. / 340 gm

**Description:** Award Winning Infant Warming Cocoon is designed for the most precious cargo. Six warming pads to cover the vital organs and two adjustable hook and loop fasteners for size adjustment and attached hood for added heat retention.

## DISASTER BLANKET

Part Number/Per Case: ED9RH/Case of 10



**Warming Temperature:**  
100°F / 38°C in 8-10 Minutes

**Maintains Heat:** 8 Hours

**Size:** 34" x 48" / 86 cm x 122 cm

**Weight:** 1.25 lbs. / .57 kg

**Description:** Emergency blanket. Used for natural disasters and mass casualty incidents. 5 year shelf life. Light and easy to store. Deploys quickly and warms immediately.

MANUFACTURED BY TECHTRADE LLC

30 Montgomery Street • Suite 690 • Jersey City, NJ 07302 • 1-866-491-TECH  
Fax: 201-706-8131 • info@ready-heat.com • www.ready-heat.com

