



THE READY-HEAT MILITARY BLANKET LINE

The Temperature Management Solution

Our proprietary panel technology is used for incremental, continuous and consistent warming. All Ready-Heat Blankets come individually sealed in an oxygen free protective poly-bag. When the bag is opened, the blanket's thermal pads react with oxygen in the atmosphere and generate heat.

Ready-Heat blankets will heat to operating temperature in 15-20 minutes and provide sufficient warmth and comfort for at least 10 hours.



Combat Casualty Care



Fixed Wing and Rotary Air Transport




Ready-Heat is part of the Hypothermia
Prevention Management Kit




4 Panel and 6 Panel Blankets

- Both blankets are US Military protocol used in the transport of wounded soldiers and the Ready-Heat 4 panel is carried by every US Army Medic.
- Due to their high tensile strength both can be used as an evacuation litter and can carry a person weighing up to 250 lbs.
- The 4 Panel is now available in a strategic military color design with blue inside and beige outside.
- Can be used in conjunction with an outer body wrap.




4 Panel

Part # S4RHMD
NSN # 6532-01-525-4062




6 Panel

Part # S6RHLG
NSN # 6532-01-525-4063



1 Panel

Part # S1RHSM
NSN # 6532-01-538-1525



12 Pad

Part # G12RH2
NSN # 6532-01-538-0759

1 Panel Blanket

- With 8 inches of attached hook & loop it can easily adjust to fit all arms and legs or be strategically placed to heat many parts of the body.
- Weighs 7 oz. (198 grams) and has a small footprint for easy addition to first aid kits.
- This blanket has been adopted by several International militaries for warming arms before administering an IV in the field.

Ready-Heat II

for High Altitude and Quick Transports

MANUFACTURED BY TECHTRADE LLC


30 Montgomery Street • Suite 690 • Jersey City, NJ 07302 • 1-866-491-TECH
Fax: 201-706-8131 • info@ready-heat.com • www.ready-heat.com



4 PANEL BLANKET

NSN #: 6532-01-525-4062

Part Number / Per Case: S4RHMD / Case of 8



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: 10 Hours

Size: 34" x 48" / 86 cm x 122 cm


Weight: 1.75 lbs. / 794 gm

Description: Torso Sized Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250 lbs./115 kg). Military Medical Protocol.

6 PANEL BLANKET

NSN #: 6532-01-525-4063

Part Number / Per Case: S6RHLG / Case of 8



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: 10 Hours

Size: 34" x 60" / 86 cm x 152 cm

Weight: 2.5 lbs. / 1.13 kg

Description: Full Body Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250lbs./115 kg). Number 1 choice of US Air Force transports.

1 PANEL BLANKET

NSN #: 6532-01-538-1525

Part Number / Per Case: S1RHSM / Case of 12



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: 10 Hours

Size: 16" x 18" / 41 cm x 46 cm


Weight: 7 oz. / 198 gm

Description: Mini Blanket. Suitable for injured extremities. Fully adjustable with 2 hook and loop closures. Perfect for First Aid and Survival Kits.

READY-HEAT II BLANKETS

NSN #: 6532-01-538-0759

Part Number / Per Case: G12RH2 / Case of 8*



Warming Temperature:
100°F / 38°C in 8-10 Minutes

Maintains Heat: 8 Hours

Size: 34" x 48" / 86 cm x 122 cm

Weight: 1.5 lbs. / 794 gm


Description: First Responder Blanket. Used to Treat Shock, Trauma and Pre-Hypothermia. Has a foam core for water resistance and is lightweight and comfortable. Quick to reach operating temperature. Works well at high altitudes.



* Also available in Torso Size
G6RH2T / Case of 12
34" x 24"/86 cm x 61 cm, 14 oz./397 gm

TEMPERATURE MANAGEMENT BLANKETS

Part Number / Per Case: SB9RH9120 / Case of 10*



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: 10 Hours

Size: 34" x 60" / 86 cm x 152 cm

Weight: 18 oz. / 510 gm

Description: Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Prevents hypothermia, increases patient comfort for higher quality of care, effective solution for infection control.



* Also available in Half Body
SB6RH9120 / Case of 20
34" x 24"/86 cm x 60 cm, 12 oz./350 gm

READY-HEAT PERSONAL USE

Part Number / Per Case: GRHV06A / Case of 12



Warming Temperature:
100°F / 38°C in 10-15 Minutes

Maintains Heat: 8 Hours


Size: Overhead Style with Adjustable Hook and Loop Straps. One Size Fits All.

Weight: 1.75 lbs. / 794 gm

Description: Designed to keep the First Responder warm during long missions. Also great for outdoor activities (Hiking, Skiing, Camping, Hunting, Fishing, Marathons and Mountain Climbing).

INFANT WARMING MATTRESS

Part Number / Per Case: S2RHIM / Case of 12



Warming Temperature:
100°F / 38°C in 15-20 Minutes

Maintains Heat: 8 Hours


Size: 11" x 13" / 28 cm x 33 cm

Weight: 7 oz. / 198 gm

Description: Infant Temperature Management. Used in conjunction with the Infant Warming Cocoon for constant warming. Great for infant transport within hospitals.

INFANT WARMING COCOON

Part Number / Per Case: GIW6C / Case of 12



Warming Temperature:
100°F / 38°C in 15 Minutes

Maintains Heat: 8 Hours


Size: 31.5" x 28" / 80 cm x 71 cm

Weight: 12 oz. / 340 gm

Description: Award Winning Infant Warming Cocoon is designed for the most precious cargo. Six warming pads to cover the vital organs and two adjustable hook and loop fasteners for size adjustment and attached hood for added heat retention.

DISASTER BLANKET

Part Number/Per Case: ED9RH/Case of 10



Warming Temperature:
100°F / 38°C in 8-10 Minutes

Maintains Heat: 8 Hours

Size: 34" x 48" / 86 cm x 122 cm

Weight: 1.25 lbs. / .57 kg

Description: Emergency blanket. Used for natural disasters and mass casualty incidents. 5 year shelf life. Light and easy to store. Deploys quickly and warms immediately.

MANUFACTURED BY TECHTRADE LLC

30 Montgomery Street • Suite 690 • Jersey City, NJ 07302 • 1-866-491-TECH
Fax: 201-706-8131 • info@ready-heat.com • www.ready-heat.com

