


## TEMPERATURE MANAGEMENT BLANKETS

### FEATURES & BENEFITS

- **Maximum Mobility:** Moves with the patient from pre-op through intra-op to post-op
- **Cost-Effective:** Only ONE blanket needed for each patient
- **Versatility:** Treats and helps prevent shock, trauma, and hypothermia in a wide variety of settings
- **Quality Materials:** Precision engineered thermal pads, medical grade bi-component non-woven material and proprietary high frequency welding
- **Air-Activated:** No electricity required
- **Warms In Minutes** to 104°F/40°C and maintains temperature for 10+ hours
- **Infection Control:** Self-heating system does not circulate particulates and is fully disposable
- **US Military Medical Protocol** includes Ready-Heat blankets as part of Tactical Combat Casualty Care since 2006
- **Five Year Shelf Life** and compact storage




#### Full Body Blanket

**Part No.:** SB9RH9120

**Open Size:** 34" x 60" (86 cm x 152 cm)

**Description:** Blue medical grade non-woven material with 9 oxygen activated engineered thermal pads.



#### Half Body Blanket

**Part No.:** SB6RH9120

**Open Size:** 24" x 34" (60 cm x 86 cm)

**Description:** Blue medical grade non-woven material with 6 oxygen activated engineered thermal pads.



#### Infant Mattress

**Part No.:** S2RHIM

**Open Size:** 11" x 13" (28 cm x 33 cm)

**Description:** Infant temperature management. Color coded blue and beige medical grade non-woven material with one oxygen activated heating panel.

MANUFACTURED BY TECHTRADE LLC


30 Montgomery Street • Suite 690 • Jersey City, NJ 07302 • 1-866-491-TECH  
Fax: 201-706-8131 • info@ready-heat.com • www.ready-heat.com



## 4 PANEL BLANKET

NSN #: 6532-01-525-4062

Part Number / Per Case: S4RHMD / Case of 8



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 34" x 48" / 86 cm x 122 cm


**Weight:** 1.75 lbs. / 794 gm

**Description:** Torso Sized Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250 lbs./115 kg). Military Medical Protocol.

## 6 PANEL BLANKET

NSN #: 6532-01-525-4063

Part Number / Per Case: S6RHLG / Case of 8



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 34" x 60" / 86 cm x 152 cm

**Weight:** 2.5 lbs. / 1.13 kg

**Description:** Full Body Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250lbs./115 kg). Number 1 choice of US Air Force transports.

## 1 PANEL BLANKET

NSN #: 6532-01-538-1525

Part Number / Per Case: S1RHSM / Case of 12



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 16" x 18" / 41 cm x 46 cm


**Weight:** 7 oz. / 198 gm

**Description:** Mini Blanket. Suitable for injured extremities. Fully adjustable with 2 hook and loop closures. Perfect for First Aid and Survival Kits.

## READY-HEAT II BLANKETS

NSN #: 6532-01-538-0759

Part Number / Per Case: G12RHZ / Case of 8\*



**Warming Temperature:**  
100°F / 38°C in 8-10 Minutes

**Maintains Heat:** 8 Hours

**Size:** 34" x 48" / 86 cm x 122 cm

**Weight:** 1.5 lbs. / 794 gm


**Description:** First Responder Blanket. Used to Treat Shock, Trauma and Pre-Hypothermia. Has a foam core for water resistance and is lightweight and comfortable. Quick to reach operating temperature. Works well at high altitudes.



\* Also available in Torso Size  
G6RH2T / Case of 12  
34" x 24"/86 cm x 61 cm, 14 oz./397 gm

## TEMPERATURE MANAGEMENT BLANKETS

Part Number / Per Case: SB9RH9120 / Case of 10\*



**Warming Temperature:**  
104°F / 40°C in 15-20 Minutes

**Maintains Heat:** 10 Hours

**Size:** 34" x 60" / 86 cm x 152 cm

**Weight:** 18 oz. / 510 gm

**Description:** Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Prevents hypothermia, increases patient comfort for higher quality of care, effective solution for infection control.



\* Also available in Half Body  
SB6RH9120 / Case of 20  
34" x 24"/86 cm x 60 cm, 12 oz./350 gm

## READY-HEAT PERSONAL USE

Part Number / Per Case: GRHV06A / Case of 12



**Warming Temperature:**  
100°F / 38°C in 10-15 Minutes

**Maintains Heat:** 8 Hours

**Size:** Overhead Style with Adjustable Hook and Loop Straps. One Size Fits All.

**Weight:** 1.75 lbs. / 794 gm

**Description:** Designed to keep the First Responder warm during long missions. Also great for outdoor activities (Hiking, Skiing, Camping, Hunting, Fishing, Marathons and Mountain Climbing).

## INFANT WARMING MATTRESS

Part Number / Per Case: S2RHIM / Case of 12



**Warming Temperature:**  
100°F / 38°C in 15-20 Minutes

**Maintains Heat:** 8 Hours


**Size:** 11" x 13" / 28 cm x 33 cm

**Weight:** 7 oz. / 198 gm

**Description:** Infant Temperature Management. Used in conjunction with the Infant Warming Cocoon for constant warming. Great for infant transport within hospitals.

## INFANT WARMING COCOON

Part Number / Per Case: GIW6C / Case of 12



**Warming Temperature:**  
100°F / 38°C in 15 Minutes

**Maintains Heat:** 8 Hours


**Size:** 31.5" x 28" / 80 cm x 71 cm

**Weight:** 12 oz. / 340 gm

**Description:** Award Winning Infant Warming Cocoon is designed for the most precious cargo. Six warming pads to cover the vital organs and two adjustable hook and loop fasteners for size adjustment and attached hood for added heat retention.

## DISASTER BLANKET

Part Number/Per Case: ED9RH/Case of 10



**Warming Temperature:**  
100°F / 38°C in 8-10 Minutes

**Maintains Heat:** 8 Hours

**Size:** 34" x 48" / 86 cm x 122 cm

**Weight:** 1.25 lbs. / .57 kg

**Description:** Emergency blanket. Used for natural disasters and mass casualty incidents. 5 year shelf life. Light and easy to store. Deploys quickly and warms immediately.

MANUFACTURED BY TECHTRADE LLC

30 Montgomery Street • Suite 690 • Jersey City, NJ 07302 • 1-866-491-TECH  
Fax: 201-706-8131 • info@ready-heat.com • www.ready-heat.com

