

Need customer support?
Scan the QR code and fill out
a support ticket to receive a solution.

Or call: (863) 414-2893

Date	Maintenance Performed	Staff Member	Notes

Suggested Maintenance and Associated Intervals:

General Inspection/operation *every 30 days (Once per month)
Drop lubrication** every three months (90 days)
Leg Bolt Tightening every six months
Kick drop down rod tightening every six months
Elevation headpiece bolt inspection every six months

*General Inspection should include:

- 1. Drop function
- 2. Condition of upholstered parts
- 3. Inspection of legs, frame and arm rests for loose attachment bolts/cracks/ damage
- 4. Operation of all features installed on the table
- **Drops should be lubricated with a single spray of WD40 on the drop plunger.