




# Measuring Chart for Mediven 550 ready to wear

medi Australia Pty Ltd

Phone: (02) 9890 8696

Fax: (02) 9890 8439

**AD** **AD Petite**



CIRCUMFERENCE IN CM											
sizes	IV (4)	V (5)	VI (6)	VII (7)	VIII (8)	IX (9)	X (10)	XI (11)	XII (12)	XIII (13)	XIV (14)
<b>cD</b>	29-33	30-34	31-35	32-36	33-37	34-38	35-39	36-40	37-41	38-42	39-43
<b>cC</b>	30-34	31-35	32-36	33-37	34-38	35-39	36-40	37-41	38-42	39-43	40-44
<b>cB1</b>	24-26	25-27	26-28	27-29	28-30	29-32	29-33	30-34	31-35	32-36	33-37
<b>cB</b>	19-20	20-21	21-22	22-23	23-24	24-25	25-26	26-27	27-28	28-29	29-30
<b>cY</b>	27-29	28-30	29-32	29-33	30-34	31-35	32-36	33-37	34-38	35-39	36-40
<b>cA</b>	18-20	19-21	20-22	21-23	22-24	23-25	24-26	25-27	26-28	27-29	28-30

**LENGTH MEASUREMENTS**

AD Regular 39-41cm

AD Petite 37-39cm

**Specials not available**

Mediven 550 Ready to wear CI.II Caramel Open Toe	
Style	Article Number
A-D	B24070
A-D Petite	B24770
A-D with narrow silicone band	B24370
A-D petite with narrow silicone band	B24C70
A-D with silicone band	B24170
A-D petite with silicone band	B24A70

**Specials not available**

Mediven 550 Ready to wear CI.II Caramel Closed Toe	
Style	Article Number
A-D	B24071
A-D Petite	B24771
A-D with narrow silicone band	B24371
A-D petite with narrow silicone band	B24C71
A-D with silicone band	B24171
A-D petite with silicone band	B24A71