


circaid® juxtafit Premium lower Leg - Size Chart

c Size	S	M	M XWide	L	L XWide	XL	XL XWide	XXL
cC	30 - 38	36 - 43	43 - 50	41 - 48	48 - 56	41 - 51	51 - 61	53 - 64
cB1	20 - 28	25 - 33	33 - 40	30 - 38	38 - 46	33 - 41	41 - 51	43 - 53
cB	15 - 23	20 - 28	20 - 28	25 - 33	25 - 33	30 - 38	30 - 38	33 - 41

Length Measurement	Short	Long
Ankle - knee crease	< 36	> 36



Measure in cm on the skin standing using a tape measure

c = Measure the circumference around the area

l = Measure the length from ankle bone to knee crease