

circaid® juxtalite®

Adjustable compression garment for treating venous diseases after initial oedema decongestion.

medi tip*:

Due to product features we recommend circaid juxtalite for:

- Women / men
- Patients who have challenges with medical compression stockings
- 24 hour compression therapy
- Mobile / immobile patients

Indications

- Chronic venous insufficiency
- Varicose veins
- Post Thrombotic Syndrome
- Venous Stasis Ulcer / Mixed Leg Ulcer
- Post Sclerotherapy
- Dependent oedema

Product features

- Easy handling – simple donning and doffing
- circaid Built-In-Pressure System (BPS): adjustable, controllable and re-adjustable compression
- Breathe-O-Prene: breathable material with high working pressure
- Latex-free, odour-inhibiting material with anti-bacterial effect
- Re-adjustability of the individual bands avoids pressure loss and thus sliding down of the garment^{1,2}
- **New:** Also available with new compressive undersock
 - Applies graduated compression to the foot and ankle
 - Graduated compression continues into the leg portion of the undersock

Benefits

- Offering an alternative therapy concept for patients who are not able to apply or tolerate ulcer compression stockings
- Faster healing rates of venous leg ulcer^{1,2}
- Promotes self-management thanks to easy handling¹


Adjustable compression ranges

20 – 30 mmHg

30 – 40 mmHg

40 – 50 mmHg

Product benefits



Built-In-Pressure



Breathe-O-Prene



Juxtaposed band



Variety



Easy handling



Inelastic

Colours

S = Standard

Beige
S

Note

* Please consider, that this recommendation does not replace the diagnosis and consultation of a medical healthcare professional.

Studies

¹ Mosti G et al. Adjustable compression wrap devices are cheaper and more effective than inelastic Bandages for venous leg ulcer healing. A Multicentric Italian Randomized Clinical Experience. *Phlebology* 2020;35(2):124-133.

² Stather PW et al. Review of adjustable velcro wrap devices for venous ulceration. In *Wound J* 2019;16(4):903-908.



circaid® juxtalite®

Your order includes

- 1 inelastic adjustable compression garment
- 1 circaid Built-In-Pressure System (BPS) guide card
- 1 pair of circaid compression anklets
- 1 circaid undersleeve lower leg




circaid juxtalite

circaid juxtalite			
Styles		AD	AD
			
		short	long
Size / Article No.	S	JU15702000	JU15002000
	M	JU15703000	JU15003000
	L	JU15704000	JU15004000
	XL	JU15705000	JU15005000
	XXL	JU15706000	JU15006000
	M extra wide	JU15713000	JU15013000
	L extra wide	JU15714000	JU15014000
	XL extra wide	JU15715000	JU15015000

Measurement chart circaid® juxtalite®

Circumferences in cm				
	S	M	M - X	L
cC	26 - 36	33 - 44	44 - 54	41 - 51
cB	19 - 22	22 - 27	22 - 27	25 - 32
	L - X	XL	XL - X	XXL
cC	51 - 61	37 - 48	48 - 61	53 - 64
cB	25 - 32	30 - 37	30 - 37	35 - 42
	Short	Long		
Sizing measurement in cm				
ℓfloor - knee crease	< 44	> 44		





circaid® juxtalite® with circaid® compressive undersock

Your order includes

- 1 inelastic adjustable compression garment
- 1 circaid Built-In-Pressure System (BPS) guide card
- 1 pair of circaid compressive undersocks moderate




circaid juxtalite with circaid compressive undersock

NEW: circaid juxtalite with compressive undersock			
Styles		AD	AD
			
		short	long
Size / Article No.	S	JU1M702000	JU1M002000
	M	JU1M703000	JU1M003000
	L	JU1M704000	JU1M004000
	XL	JU1M705000	JU1M005000
	XXL	JU1M706000	JU1M006000
	M extra wide	JU1M713000	JU1M013000
	L extra wide	JU1M714000	JU1M014000
	XL extra wide	JU1M715000	JU1M015000

Measurement chart circaid® juxtalite®

Circumferences in cm				
	S	M	M - X	L
cC	26 - 36	33 - 44	44 - 54	41 - 51
cB	19 - 22	22 - 27	22 - 27	25 - 32
	L - X	XL	XL - X	XXL
cC	51 - 61	37 - 48	48 - 61	53 - 64
cB	25 - 32	30 - 37	30 - 37	35 - 42
	Short		Long	
Sizing measurement in cm				
f _{floor - knee crease}	< 44		> 44	



circaid[®] foot options

The following foot options are recommended to wear in combination with circaid juxtalite.

Note

It is recommended to always wear a circaid foot option in combination with a circaid leg garment. The indications of the circaid foot options are identical to the indications of the corresponding circaid leg garment. Please order separately.

circaid[®] juxtalite[®] ankle foot wrap

Adjustable compression device with closed heel for venous diseases.

Product details


- Provides compression to foot and ankle
- Can be trimmed in length to provide perfect fit
- Easy handling
- Available in standard sizes only

Measurement chart

Circumferences in cm		
	S	M
cA1	19 – 26	24 – 31
Article No.	JU1W002000	JU1W003000



circaid juxtalite ankle foot wrap



circaid[®] compression anklet

Standard compression option for foot and ankle.

Product details


- Provides compression to the foot and ankle
- Can be combined with other circaid foot options
- Available in standard and large

Measurement chart

Circumferences in cm		
	Standard	Large
cB	19 – 27	25 – 50
cA1	19 – 36	22 – 41
Article No.	JU51001000	JU51002000



circaid compression anklet



circaid® power added compression band (pac band)

Additional compression for the arch of the foot and ankle.

Product details


- Additional compression to the foot with mild to severe tendency for oedema
- Targets the arch of the foot and around the malleoli
- Can be trimmed to fit almost every foot shape
- Easy handling
- Lightweight and breathable material
- Available in standard and large
- Bilateral use (left or right foot)

Measurement chart

Circumferences in cm			
	Standard	Large	Extra Large
cB	19–27	25–50	45–56
cA1	19–36	22–41	22–41
Article No.	JU5P001000	JU5P002000	JU5P003000



circaid power added compression band (pac band)



circaid® compressive undersock

Easy-to-apply liner with targeted compression for foot and ankle.

Product details

- Applies graduated compression to the foot and ankle
- Graduated compression continues into the leg portion of the undersock
- Soft material for pleasant wearing comfort
- Available in 15–25 mmHg (moderate compression) and 25–35 mmHg (firm compression)
- Available in standard and large

Measurement chart

Circumferences in cm		
	Standard	Large
cA1	19–32	24–41
cB	19–27	25–42
cC	26–54	41–64
Article No.		
Moderate	JU5M001000	JU5M002000
Firm	JU5F001000	JU5F001000



circaid compressive undersock

