




Measuring Chart for Mediven 550 ready to wear

medi Australia Pty Ltd

Phone: (02) 9890 8696

Fax: (02) 9890 8439

AD **AD Petite**



CIRCUMFERENCE IN CM											
sizes	IV (4)	V (5)	VI (6)	VII (7)	VIII (8)	IX (9)	X (10)	XI (11)	XII (12)	XIII (13)	XIV (14)
cD	29-33	30-34	31-35	32-36	33-37	34-38	35-39	36-40	37-41	38-42	39-43
cC	30-34	31-35	32-36	33-37	34-38	35-39	36-40	37-41	38-42	39-43	40-44
cB1	24-26	25-27	26-28	27-29	28-30	29-32	29-33	30-34	31-35	32-36	33-37
cB	19-20	20-21	21-22	22-23	23-24	24-25	25-26	26-27	27-28	28-29	29-30
cY	27-29	28-30	29-32	29-33	30-34	31-35	32-36	33-37	34-38	35-39	36-40
cA	18-20	19-21	20-22	21-23	22-24	23-25	24-26	25-27	26-28	27-29	28-30

LENGTH MEASUREMENTS

AD Regular 39-41cm

AD Petite 37-39cm

Specials not available

Mediven 550 Ready to wear CI.II Caramel Open Toe	
Style	Article Number
A-D	B24070
A-D Petite	B24770
A-D with narrow silicone band	B24370
A-D petite with narrow silicone band	B24C70
A-D with silicone band	B24170
A-D petite with silicone band	B24A70


Specials not available

Mediven 550 Ready to wear CI.II Caramel Closed Toe	
Style	Article Number
A-D	B24071
A-D Petite	B24771
A-D with narrow silicone band	B24371
A-D petite with narrow silicone band	B24C71
A-D with silicone band	B24171
A-D petite with silicone band	B24A71

Measuring Chart for Mediven Mondi ready to wear



AD AD Petite AG PrB AG Petite PrB



CIRCUMFERENCE IN CM							
sizes	I	II	III	IV	V	VI	VII
cG	44-43	46-57	50-61	53-65	56-70	60-74	64-80
cF	40-47	42-51	45-55	48-58	51-62	54-67	58-72
cE	30-37	33-40	36-42	38-44	40-46	42-48	44-50
cD	28-33	30-35	32-37	34-39	36-41	38-43	40-45
cC	29-34	31-36	33-38	35-40	37-42	39-44	41-46
cB1	24-26	26-28	28-31	30-33	32-35	34-37	36-39
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cY	28-30	29-32	31-34	33-36	34-38	36-40	37-41
cA	18-21	20-23	22-25	24-27	26-29	28-31	30-33

LENGTH MEASUREMENTS
 AD Regular 38-43cm AG Regular 72-83cm
 AD Petite 33-37cm AG Petite 62-71cm

Specials not available

Mediven Mondi Ready to Wear Cl.II Caramel Open Toe	
Style	Article Number
A-D	C24070
A-D Petite	C24770
A-G Petite with silicone topband	C25170
A-G with silicone topband	C26170
A-D with narrow silicone band	C24170
A-D petite with narrow silicone band	C24A70
A-D with silicone band	C24370
A-D petite with silicone band	C24C70

Specials not available

Mediven Mondi Ready to Wear Cl.II Caramel Closed Toe	
Style	Article Number
A-D	C24071
A-D Petite	C24771
A-D with narrow silicone band	C24171
A-D petite with narrow silicone band	C24A71
A-D with silicone band	C24371
A-D petite with silicone band	C24C71