

Mediven Espirit Ready to wear Armsleeves & Gloves

AG Combined Armsleeve	AG Combined Armsleeve with Topband
--------------------------	---------------------------------------




Standard
 CCL1: J15070x
 CCL2: J25070x
Long
 CCL1: J16070x
 CCL2: J26070x



Standard
 CCL1: J15170x
 CCL2: 555 J25170x
Long
 CCL1: J16170x
 CCL2: 556 J26170x

CG Armsleeve	CG Armsleeve with Topband	AC1 glove with compressive fingers
-----------------	------------------------------	---------------------------------------



Standard
 CCL1: J1H070x
 CCL2: J2H070x
Long
 CCL1: J1G070x
 CCL2: J2G070x



Standard
 CCL1: J1H170x
 CCL2: 550 J2H170x
Long
 CCL1: J1G170x
 CCL2: 551 J2G170x



CCL1: J11071x
 CCL2: 552 J21071x


x=size (1-5)

RAL CCL1: (15-21 mmHg)

Caramel only

RAL CCL2: (23-32 mmHg)


Mediven Espirit Ready to wear Armsleeves & Gloves sizing chart



Mediven Espirit Glove sizing chart

Sizes	I	II	III	IV	V
cA	17.0-19.0	18.0-20.0	19.0-22.0	21.0-24.0	23.0-27.0
cB	18.0-20.0	19.0-21.0	20.0-23.0	22.0-25.0	24.0-28.0
cC	15.0-17.0	17.0-19.0	19.0-21.0	21.0-23.0	23.0-25.0
cC1	16.0-19.0	18.0-21.0	20.0-23.0	22.0-25.0	24.0-27.0

Circumference/Fingers (cm)




Finger		I	II	III	IV	V
Thumb	top cZ	5.5	6.5	7.5	8.5	9.5
	bottom cX	6.5	7.5	8.5	9.5	10.5
Index	top cZ	4.5	5.5	6.5	7.5	8.5
	bottom cX	5.5	6.5	7.5	8.5	9.5
Middle	top cZ	4.5	5.5	6.5	7.5	8.5
	bottom cX	5.5	6.5	7.5	8.5	9.5
Ring	top cZ	4	5	6	7	8
	bottom cX	5	6	7	8	9
Little	top cZ	3.5	4.5	5.5	6.5	7.5
	bottom cX	4.5	5.5	6.5	7.5	8.5

Length/Fingers (cm)

Thumb	3.5-4.5
Index	4.5-5.5
Middle	5.5-6.5
Ring	4.5-5.5
Little	3.5-4.5

Length/Handpart (cm) all sizes

AB	3.5-4.5
AC	10.0-12.0
AC1	15.0-16.0



Mediven Espirit Armsleeve sizing chart

Sizes	I	II	III	IV	V
cA	17.0-19.0	18.0-20.0	19.0-22.0	21.0-24.0	23.0-27.0
cC	15.0-17.0	17.0-19.0	19.0-21.0	21.0-23.0	23.0-25.0
cE	25.0-27.5	27.0-29.5	29.0-32.0	31.0-34.0	33.0-37.0
cG	28.0-32.0	30.0-34.0	32.0-37.0	35.0-40.0	38.0-44.0

Length cG: Standard 40-44cm Long 44-48cm

