


mediven® comfort Below Knee



Measure on the skin while standing
Use a tape measure

Length - Measure from floor up to:
A to D (Below Knee)

Size - Measure around the area:
B (Ankle)
C (Calf)
D (Below Knee)

SIZE	1	2	3	4	5	6	7
cD	27 - 33	29 - 36	32 - 39	34 - 42	36 - 45	38 - 48	40 - 50
cC	28 - 34	30 - 37	33 - 40	35 - 43	37 - 46	39 - 49	41 - 51
cB	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30	30 - 32

c = Circumference in cm

	Measure	petite	regular
Length in cm	ID (AD)	34 - 38	39 - 44

l = length in cm