

Product benefits



Clima Comfort



Clima Fresh



Perfect Fit



cf. Ben Ammar, F. (Dermiscan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

mediven for men

Styles

AD

AD pte.



calf

calf petite

Class	Article No.	3A40x11-7	3A47x11-7
Class 1	Article No.	3A40x11-7	3A47x11-7
Class 2	Article No.	3B40x11-7	3B47x11-7

See the last page for the article numbers

Standard styles

Standard	Size I – VII, Length standard and petite
Special models	not available
Border	AD Tricot border
Foot part	closed toe shoe size 39 – 46

Material composition

	Elastane (EL)	Polyamide (PA)
CCL 1	23%	77%
CCL 2	29%	71%

Material

Light

Compression class

CCL 1, 2

Standard colours

S = Standard

Black
S

Anthracite
S

Navy
S

Brown
S



medi tip:



medi Butler donning made easier¹ medi skin and textile care offer support

More information on the medi accessories from page 112

¹ In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.