

mediven® forte Below Knee - Standard Size Chart

SIZE	1	2	3	4	5	6	7
cD	28 - 33	30 - 35	32 - 37	34 - 39	36 - 41	38 - 43	40 - 45
cC	29 - 34	31 - 36	33 - 38	35 - 40	37 - 42	39 - 44	41 - 46
cB	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30	30 - 32

c = Circumference in cm l = length in cm

Measure Length	petite	regular
Floor to D	33 - 37	38 - 43

Measure on the skin while standing using a tape measure

Size - Measure around the area B (Ankle) C (Calf) D (Below Knee)

Length - Measure from floor to D (Below Knee)

