


# mediven® forte Thigh Stocking + Silicone Top Band Size Chart



| SIZE | 1       | 2       | 3       | 4       | 5       | 6       | 7       |
|------|---------|---------|---------|---------|---------|---------|---------|
| cG   | 44 - 47 | 46 - 50 | 50 - 54 | 53 - 58 | 56 - 61 | 60 - 65 | 64 - 70 |
| cD   | 28 - 33 | 30 - 35 | 32 - 37 | 34 - 39 | 36 - 41 | 38 - 43 | 40 - 45 |
| cC   | 29 - 34 | 31 - 36 | 33 - 38 | 35 - 40 | 37 - 42 | 39 - 44 | 41 - 46 |
| cB   | 18 - 20 | 20 - 22 | 22 - 24 | 24 - 26 | 26 - 28 | 28 - 30 | 30 - 32 |

| Length | Petite  | Regular |
|--------|---------|---------|
| IG     | 62 - 71 | 72 - 83 |

Measure on the skin while standing - Use a tape measure

c = Circumference in cm  
l = Length in cm

Size - Measure around the area  
B (Ankle)  
C (Calf)  
D (Knee)  
G (Thigh)

Length Measure from floor up to G (Thigh)