

mediven® Measurement charts

Measurement chart

mediven® forte

Standard fit without seam		I	II	III	IV	V	VI	VII	
AG with topband, normal	cG	44-47	46-50	50-54	53-58	56-61	60-65	64-70	
AG with topband, wide	cG	48-53	51-57	55-61	59-65	62-70	66-74	71-80	
	cG	44-53	46-57	50-61	53-65	56-70	60-74	64-80	
	cF	40-47	42-51	45-55	48-58	51-62	54-67	58-72	
	cE	30-37	33-40	36-42	38-44	40-46	42-48	44-50	
	cD	28-33	30-35	32-37	34-39	36-41	38-43	40-45	
	cC	29-34	31-36	33-38	35-40	37-42	39-44	41-46	
	cB1	24-26	26-28	28-31	30-33	32-35	34-37	36-39	
	cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32	
	cY	28-30	29-32	31-34	33-36	34-38	36-40	37-41	
	cA	18-21	20-23	22-25	24-27	26-29	28-31	30-33	
		Length measurements							
cG regular		72-83 cm		Foot					
cG petite		62-71 cm		2 foot lengths, closed toe					
dD regular		38-43 cm		Small foot up to 23 cm ≙ Shoe size 36-38					
dD petite	33-37 cm		Large foot from 23 cm and up ≙ Shoe size 38-43						

Measurement chart

mediven active® · mediven® for men mediven® angio

Standard fit without seam		I	II	III	IV	V	VI	VII	
	cD	27-33	29-36	32-39	34-42	36-45	38-48	40-50	
	cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51	
	cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41	
	cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32	
	cY	26-31	28-33	29-35	31-37	32-38	33-40	34-42	
	cA	17-22	19-24	21-26	23-29	25-32	27-34	28-35	
	Length measurements								
	mediven active			Foot					
	dD regular	43-49 cm		mediven active		closed toe		shoe size 38-48	
	dD petite	37-42 cm		mediven for men		closed toe		shoe size 39-46	
mediven for men			mediven angio		mediven angio		2 foot lengths, closed toe		
dD regular	39-44 cm				Small foot		up to 25 cm ≙ Shoe size 36-40		
dD petite	34-38 cm				Large foot		from 25 cm ≙ Shoe size 40-46		