


mediven harmony® Size Chart



Standard Size

Length C to G


Standard 41 - 47cm

All measurements are in cm with Tape Measure

1. In a relaxed position without tension place your arm on a table to measure the circumference of: **C, C1, D, E**
2. Measurement **cC1** is 5cm above **cC**
3. With a slightly bent arm measure **E**
4. While standing measure **F & G** with arm beside your body
5. **Length** is taken measuring the inside of a straightened arm

SIZE CM	1	2	3	4	5	6	7
cG	25.5 - 28	28 - 30.5	30 - 33	32 - 35	33.5 - 36.5	36 - 39.5	38 - 41.5
cF	23.5 - 25.5	25 - 27.5	27.5 - 30	29.5 - 32	31 - 34	33 - 36	35 - 38.5
cE	23 - 25	24.5 - 26.5	26 - 28.5	27.5 - 30	29 - 31.5	30.5 - 33.5	32.5 - 35.5
cD	22 - 24	23 - 25	24 - 26	25 - 27.5	27 - 29.5	28.5 - 31	30.5 - 33.5
cC1	15 - 16.5	16.5 - 18	18 - 19.5	19.5 - 21	21 - 23	23 - 25	25 - 27.5
cC	13.5 -14.5	14.5 -15.5	16 - 17.5	17.5 - 19	19 - 20.5	20.5 -22.5	22.5 -24.5
cA	15 - 18	16 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 27

mediven harmony® Size Chart



Extra Wide Size

Length C to G

Standard 41 - 47cm

All measurements are in cm with Tape Measure

1. In a relaxed position without tension place your arm on a table to measure the circumference of: **C, C1, D, E**
2. Measurement **cC1** is 5cm above **cC**
3. With a slightly bent arm measure **E**
4. While standing measure **F & G** with arm beside your body
5. **Length** is taken measuring the inside of a straightened arm

SIZE CM	1	2	3	4	5	6	7
cG	28.5 - 31	30.5 - 33.5	33 - 36	35 - 38.5	37 - 40.5	39.5 - 43	41.5- 45.5
cF	25.5 -28	27.5 - 30	30 - 33	32 - 35	34 - 37	36 - 39.5	38.5 - 42
cE	24 - 26	26.5 - 29	28.5 - 31	30 - 33	31.5 - 34.5	33.5 - 36.5	35.5 - 39
cD	24 - 26	25 - 27.5	26 - 28.5	28 - 30.5	29.5 - 32	31 - 34	33.5 - 36.5
cC1	15 - 16.5	16.5 - 18	18 - 19.5	19.5 - 21	21 - 23	23 - 25	25 - 27.5
cC	13.5 -14.5	14.5 -15.5	16 - 17.5	17.5 - 19	19 - 20.5	20.5 -22.5	22.5 -24.5
cA	15 - 18	16 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 27