



01928 262128
info@dentaledge.co.uk
www.dentaledge.co.uk

Hassle-free Engineering

Is your equipment in need of service or repair?
Did you know Dental Edge UK will let you know
when service and inspection work is due?

**We make the process easy and
incredibly convenient.**



A COMPLETE SERVICE

From Repair
to Validation.



SOUTH WEST PRESENCE

Less than 1.5 hrs
travel to most
customers.



MULTI SKILLED ENGINEERS

Fully stocked and
manufacturer
trained.



HUGE SELECTION OF PARTS

Stocked in our
warehouse and
vans.

Engineering Services.....

Dental Edge UK offer installation, service, repair and technical support for all major manufacturers. For sales or aftersales get in touch with a member of our friendly team who will be happy to help.

SERVICING



Autoclaves



Compressor's



Digital Imaging



Reverse Osmosis
Systems



Suction System's



Treatment
Centre's



Ultrasonic Bath's



Washer
Disinfectors

VALIDATION



Autoclaves



Ultrasonic Bath's



Washer
Disinfectors

INSPECTION



Autoclaves



Compressors



X-Ray's

REPAIR

We understand the urgency and importance of getting your equipment back up and running. With our experience and dedication we will ensure little or no downtime, no loss of business, and a fast and efficient service. Our engineers will give free training and advice on how to best maintain your equipment to reduce costly repairs.

Health and Social Care Act 2008 (Regulated Activities) Regulations 2014: Regulation 15

"All equipment must be used, stored and maintained in line with manufacturers' instructions."

We can help with the key legislation that dental organisations are obliged to comply with.

See next page.

Equipment Inspection

Pressure Vessel Inspection

Written Scheme of Examination

- Compressor Pressure Vessel Inspection (PVI)
- Autoclave Pressure Vessel Inspection (PVI) (N,B and S type)

Pressure Systems Safety Regulations 2000 (PSSR 2000)

Provision and Use of Work Equipment Regulations (PUWER) 1992.

It is a mandatory requirement that compressors, autoclaves and their associated pipework and protective devices are regularly inspected as specified in the written scheme of examination and by a competent person. Our engineers are qualified to complete a full Written Scheme of Examination in accordance with the **Pressure Systems Safety Regulations 2000**. This in turn will form the basis of compliance to the **Provision and Use of Work Equipment Regulations (PUWER) 1992**.

Without a valid inspection you may also be invalidating your practice insurance. We will issue a reminder well in advance to ensure you remain compliant.

X-Ray Inspection

- Critical Examination (on installation)
- Routine Inspection (every 3 years)
- AXTS (annual X-ray test service)

The Ionising Radiation Regulations 2017 (IRR17)

Ionising Radiation (Medical Exposure) Regulations 2017 (IRMER17)

Referring to **Dental Guidance Notes 2020, Section 5.2.3** you must have a routine inspection of your x-ray equipment carried out at least every 3 years however, this will not adequately cover the health and safety regulations that you ensure your X-ray equipment is safe for use. Therefore you should consider an annual electrical and mechanical safety check (AXTS) to ensure the continued safe use of your equipment.

Dental Edge UK Ltd have worked closely with an RPA and MPE to comply with **IRR17** to carry out inspections from the point of installation to annual service checks to ensure you comply in the safe use of your X-ray equipment. We will issue a reminder well in advance to ensure you remain compliant.

FOR FURTHER DETAILS ON ANY SERVICE WE OFFER PLEASE VISIT OUR WEBSITE:
WWW.DENTALEDGEUK.CO.UK OR GIVE US A CALL: 01928 262128

Service comes as Standard

We strive to provide the quickest and most hassle-free engineering service available.



Why Settle for less?

We are passionate and determined to go above and beyond with service. From the moment we first talk, through to completion of work we operate on the basis of making it easy so you can do what you do best – looking after your patients.

www.dentaledgeuk.co.uk

Shop our full range of self fit spares and equipment with useful diagrams and compatibility information.



Scan with camera

Recycle Centre

Quality pre owned products and warehouse clearance stock at great prices.



Scan with camera

