

Matheson Dexterity and Range of Motion Portable Panel

Evaluate fine finger dexterity, medium dexterity, gross grasp, hand-eye coordination, sustained sitting, and sustained standing with one portable system.

This system is a criterion referenced test. Method-Time Measurement (MTM) time standards have been established for the 5 protocols.

Protocols:

- Fine Finger Placing and Displacing Protocol
- Small Screw Assembly and Disassembly
- Medium Screw Assembly and Disassembly
- Peg Place and Displace
- Peg Turn and Place

