

# Protective Face Shield

Patent Pending

**Product Features:**

- Optimal fit—can be comfortably worn for many hours
- Allows visibility for facial expressions to maintain a positive learning environment
- Redirects bi-directional airflow
- Can easily adjust up for active engagement and down for idle circumstance
- Breakaway safety feature—allows for quick attaching and detaching from neck
- Ultra lightweight
- Minimal noticeability—forms to the user’s neck and chest



FDA Registration Number 3016840036

SKU# 300-024



**Shield Bracket Dimensions:**  
3.68" x 1.37" x 2.86"



**Neck Mount Dimensions:**  
4.18" x 4.2" x 1.4"



**Assembled Shield:**  
Active State



**Assembled Shield:**  
Idle State

**Materials Used:**

Frame: Polypropylene (Black)

Visor: PETG (Optically Clear)

**Customization:**

Ask about customizing with your organizations primary color and logo!

**Cleaning Instructions**

Wash all components with soap and water


Frame—Dishwasher safe

Note: All components are shipped unsterilized

**Please Recycle**



FRAME



**For more information:**  
[covidsupport@zverse.com](mailto:covidsupport@zverse.com)  
[www.zverse.com](http://www.zverse.com)

**Order online at:**  
[BeTheShield.info](http://BeTheShield.info)

    @zverse3d  
Please share your ZShields in action using #BeTheShield

Personal Protective Equipment EUA

On the basis of this determination, the Secretary of Health & Human Services declared that circumstances exist justifying the authorization of emergency use of personal respiratory protective devices during the COVID-19 outbreak pursuant to section 564 of the FD&C Act, subject to the terms of any authorization issued under that section.

If you have any questions, please email: [CDRH-NonDiagnosticEUA-Templates@fda.hhs.gov](mailto:CDRH-NonDiagnosticEUA-Templates@fda.hhs.gov)

ZVERSE INC. ALL RIGHTS RESERVED 2020

