









VFT-28

VACCINE FRIDGE
DIGITAL THERMOMETER

SPECIFICATIONS

FEATURES



DISPLAY

Top - Current temperature

Left - Min temperature

Right - Max temperature



DIMENSIONS

Width: 98mm

Depth: 22mm

Height: 110mm



FEATURES

- Exceptionally easy to view and record hi/lo readings
- Hi, lo and current temperatures always on display
- One button press to reset hi/lo readings
- Audible alarm when outside +2°C to +8°C range
- Alarm has one press deactivation
- Perfect temeprature monitor for vaccine fridge or backup cooler
- No setup required
- Includes 2 sensors to measure air or vaccine temperature
- Designed and manufactured specifically for Australian and New Zealand vaccine fridges
- Fold out stand or hang from a hook
- 2x AAA batteries included

Three Year Warranty

Office 1800 664 664

Service 1800 630 360