

# LKexv 2600

## MediLine



<b>Total gross capacity</b>	240 litre
<b>Ext. dimensions</b> in mm (w/d/h)	600/610/1250
<b>Int. dimensions</b> in mm (w/d/h)	470/440/1062
<b>Energy consumption</b> in 24 h	0,786 kWh
Climate rating <sup>1</sup>	SN-T
Cooling system	forced-air
Defrosting	automatic
Temperature range	+1°C to +15°C
Gastronorm	–
Housing material/colour	steel/white
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	analogue controls
Temperature display	external digital
Interior light	–
Fault: warning signal	–
Interface	–
volt-free alarm contact	no
Shelves	5
Adjustable Shelving	4
Number of drawers	–
Number of baskets	–
Maximum number of baskets	–
Shelf material	glass
Shelf loading	40 kg
Height of adjustable legs	-
Castors/type of castors	–
Door handle	Ergonomic door handle
Foot pedal opener	–
Self-closing door	–
Insulation	–
Defrost water drain	–
Refrigerant	R 600a
Noise level	48 dB(A)
Gross/net weight	53/49 kg
EAN-Nr.	9005382164891
<b>Connection Rating</b>	1.0 A / 120 W

<sup>1</sup> Climate Ratings: SN +10°C to +32°C N +16°C to +32°C ST +16°C to +38°C T +16°C to +43°C