



Features:

- Fan Forced System
- TGA Approved
- Voltage Free Connection
- Door Lock
- Door Ajar Alarms
- Interior LED Lighting
- Automatic Condensate water Disposal
- Fan Door Switch
- Data Logger Included
- Battery Backup for Controller
- Digital Display with Built In Alarms
- Hi / Lo Read out

Temperature

2°C to 8°C




Capacity

ML29	29L
ML59	59L
ML125	125L
ML300	300L
ML400	428L



Dimensions



ML29	445w x 485d x 535h
ML59	445w x 485d x 725h
ML125	595w x 595d x 820h
ML300	595w x 650d x 1500h
ML400	595w x 695d x 1845h

Specifications

*Refer to warranty sheet

	ML29	ML59	ML125	ML300	ML400
Capacity	29 L	59 L	125 L	300 L	428 L
External Dimensions (including castors)	445w x 485d x 535h	445w x 485d x 725h	595w x 595d x 820h	595w x 650d x 1500h	595w x 695d x 1845h
Internal Dimensions	360w x 400d x 400h	360w x 200d x 590h	510w x 410d 680h	510w x 465d x 1360h	510w x 485d x 1725h
Weight	25 Kg	28 Kg	42 Kg	72 Kg	84 Kg
Voltage	240V 50Hz	240V 50Hz	240V 50Hz	240V 50Hz	240V 50Hz
Watts	100w / 0.8A	100w / 0.8A	135w / 1.05A	240w / 1.5A	280w / 1.7A
Refrigerant	R134a	R134a	R134a	R134a	R134a
Warranty	12 Months*	12 Months*	12 Months*	12 Months*	12 Months*
Shelves Size	348x192	(2) 348x328 & (1) 348x192	(3) 495x405 & (1) 495x260	495x450	495x450
Number of Shelves	2 Adjustable	3 Adjustable	3 Adjustable	5 Adjustable	5 Adjustable
Digital Display	✓	✓	✓	✓	✓
Castors	Optional	Optional	Optional	Optional	Optional
Interior/Exterior Cabinet Finish	White Enamel	White Enamel	White Enamel	White Enamel	White Enamel
Insulation	44mm	44mm	44mm	44mm	44mm
High/Low Temperature Alarm	✓	✓	✓	✓	✓
Fan Forced Refrigeration System	✓	✓	✓	✓	✓
Door Lock	✓	✓	✓	✓	✓
Solid Door / Glass Door	✓	✓	✓	✓	✓
Independent data logger	✓	✓	✓	✓	✓
Optional Equipment					
Door ajar alarm	✓	✓	✓	✓	✓
Port entry	✓	✓	✓	✓	✓
Remote alarm	✓	✓	✓	✓	✓
Temp alert	✓	✓	✓	✓	✓
Chart recorder	✓	✓	✓	✓	✓
Battery backup for controller	✓	✓	✓	✓	✓
Castors	✓	✓	✓	✓	✓