



Features

- Top mounted system
- TGA Approved
- Five (5) rows of PVC coated shelves behind each door
- Interior Lighting LED
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms
- Hi / Lo Read out

Temperature




10°C to 40°C

Capacity



NHRiT400 350L

NHRiT600 590L



Dimensions

NHRiT400 600w x 620d x 1915h

NHRiT600 770w x 750d x 1965h



Optional equipment

Temperature Data Logger



Keypad

Model	NHRiT400	NHRiT600
Capacity	350 L	590 L
External Dimensions (including castors)	600w x 620d x 1915h	770w x 750d x 1965h
Internal Dimensions	500w x 470d x 1490h	670w x 600d x 1530h
Weight	95 Kg	115 Kg
Voltage	240V 50Hz	240V 50Hz
Watts	670W	850W
Power Requirement	10amp GPO	10amp GPO
Refrigerant	R134a	R134a
Warranty	12 Months *	12 Months *
Shelves size	480 x 450	650 x 580
Number of Shelves	5 Adjustable	5 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50mm	50mm
High/Low Temperature Alarm	✓	✓
Fan Forced Refrigeration System	✓	✓
Door Lock	✓	✓
Glass Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temperature alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

**Mediline reserves the right to alter the specifications without prior notice. **