

NULINE
REFRIGERATION

Win a Free Medical Fridge or Freezer*



We are searching for the oldest working Nuline Refrigeration's medical fridge or freezer as part of our 40th Anniversary celebration.

Win a new Nuline Refrigeration's medical fridge most like your existing one.

Proudly Manufacturing in Australia and a locally owned family business.

Simply fill in the form on our website

Submissions close June 30th 2025, and the winning entrant will be notified via email and phone on July 14th 2025.

* Terms and Conditions apply



LABORATORY INCUBATOR

NULINE
REFRIGERATION

FEATURES

- Fan Forced System
- Overheating Alarm
- Safety Cutoff at 70°C
- Double Glazed Doors
- Stainless int./ext. with Chrome Plated Shelves
- Microprocessor Control
- Digital Display with High Alarm
- ARTG Number 254719



	Ambient +5°C to 60°C
	80 L 247 L 611 L 1000 L
	500 x 600 x 625 700 x 700 x 972 750 x 700 x 1880 1200 x 700 x 1880



LABORATORY INCUBATOR



Model	Li100	Li260	Li650	Li1100
Gross Capacity	80 L	247 L	611 L	1000 L
External Dimensions in mm (w/d/h including castors)	500 x 600 x 625	700 x 700 x 972	750 x 700 x 1880	1200 x 700 x 1880
Internal Dimensions in mm (w/d/h)	390 x 480 x 370	590 x 685 x 600	630 x 610 x 1590	1100 x 550 x 1555
Packed Dimensions in mm (w/d/h including pallet)	575 x 670 x 765	800 x 800 x 1070	850 x 800 x 2030	1300 x 800 x 2030
Gross Weight - kg (packed weight)	63	109	180	260
Net Weight - kg	48	89	160	220
Voltage	240V			
Temperature Fluctuation at 30°C	±0.1°C	±0.5°C	±0.9°C	±0.9°C
Temperature Range	Ambient +5°C to 60°C			
Ambient Temperature Range	10°C to 32°C			
Exterior Cabinet Finish	Stainless Steel			
Inner Liner Material	Stainless Steel			
Door	Glass			
Door Hinge	RH			LH / RH
Door Handle	Full Length			
Type of Controller	Microprocessor			
Temperature Display	Digital			
BMS Connectivity	Optional			
(BMS Connectivity for Remote Alarms only)				
(BMS Connectivity for 4-20mA Transmitter)				
Port Entry	Optional			
Power Requirement	10amp GPO			
Castors	Fixed Legs			Yes



Nuline Refrigeration reserves the right to alter the specifications without prior notice.

AUSTRALIAN MADE

The logo for Nuline Refrigeration is an oval shape with a blue gradient. The word "NULINE" is written in large, white, bold, sans-serif capital letters. Below it, the word "REFRIGERATION" is written in smaller, white, sans-serif capital letters.

NULINE
REFRIGERATION

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to -40°C), blanket warmer, mortuary refrigerators and refrigerated incubators.

