

# NDT 700 REFRIGERATOR AND FREEZER

**NULINE**  
REFRIGERATION

The NDT Spark Free contains a dual refrigeration system: a unit for refrigeration and another for freezer.

This locally manufactured cabinet offers quality 304 grade stainless steel interior and 430 grade stainless steel exterior as standard with the option of 316.

## FEATURES

- Top mounted system
- Heated cabinet facia to prevent condensation
- Auto defrost in 2°C Section
- Manual defrost in -20°C section
- Stainless steel interior / exterior
- Digital display with built in alarms
- Hi/Lo Readout



Refrigerator: 2°C to 8°C  
Freezer: -18°C to -22°C



Refrigerator: 386 L  
Freezer: 212 L



870w x 750d x 2120h



# NDT 700 REFRIGERATOR AND FREEZER



Model	Refrigerator	Freezer
Gross Capacity	386 L	212 L
External Dimensions in mm (w/d/h including castors)	870 x 750 x 2120	
Internal Dimensions in mm (w/d/h)	770 x 725 x 885	685 x 535 x 580
Packed Dimensions in mm (w/d/h including pallet)	920 x 800 x 2250	
Gross Weight - kg (packed weight)	132	
Net Weight - kg	125	
Voltage	240V	
Temperature Range	2°C to 8°C	-18°C to -22°C
Ambient Temperature Range	10°C to 32°C	
Refrigerant	R134a	R404a
Cooling System / Defrost	Static	
Exterior Cabinet Finish	Stainless Steel	
Inner Liner Material	Stainless Steel	
Door	Solid	
Door Hinge	LH	
Door Handle	Lockable	
Type of Controller	Digital	
Temperature Display	External	
Data Logger	Optional	
Door Ajar Alarms	Optional	
High and Low Readout	Visual / Audible	
BMS Connectivity (BMS Connectivity for Remote Alarms only)	Optional	
(BMS Connectivity for 4-20mA Transmitter)	Optional	
Port Entry	Optional	
Battery Back-up for Controller	Optional	
Power Requirement	10amp GPO	
Watts	960W	540W
Castors	Yes	
Adjustable Shelves	3	2
Usable Shelf Area in mm(w/d)	530 x 652	675 x 520



Nuline Refrigeration reserves the right to alter the specifications without prior notice.

# NDT 700 REFRIGERATOR AND FREEZER

**NULINE**  
REFRIGERATION

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to  $-40^{\circ}\text{C}$ ), blanket warmer, mortuary refrigerators and refrigerated incubators.

