




# NLMi REFRIGERATOR INCUBATOR

**NULINE**  
REFRIGERATION

## FEATURES

- Top Mounted System
- Five (5) rows of PVC Coated Shelves behind each door
- Interior Lighting LED
- Castors
- Automatic Water Disposal
- Digital Display with Built-in Alarms
- Hi/Lo Readout



	2°C to 40°C
	400 L 700 L 1000 L
	600w x 620d x 1915h 870w x 700d x 2100h 1200w x 700d x 2100h



# NLMi REFRIGERATOR INCUBATOR



Model	NLMi400/1	NLMi700/1	NLMi1000/2
Gross Capacity	400 L	700 L	1000 L
External Dimensions in mm (w/d/h including castors)	600 x 620 x 1915	870 x 700 x 2100	1200 x 700 x 2100
Internal Dimensions in mm (w/d/h)	500 x 470 x 1435	770 x 550 x 1555	1100 x 550 x 1555
Packed Dimensions in mm (w/d/h including pallet)	700 x 720 x 2065	970 x 800 x 2250	1300 x 800 x 2250
Gross Weight - kg (packed weight)	102	112	137
Net Weight - kg	95	95	125
Voltage	240V		
Temperature Range	2°C to 40°C		
Ambient Temperature Range	10°C to 32°C		
Temperature Fluctuation	±1.5°C		
Refrigerant	R134a		
Cooling System / Defrost	Fan Forced / Automatic		
Exterior Cabinet Finish	White Colourbond		
Inner Liner Material	White Colourbond		
Door	Solid / Glass E-Glass		
Door Hinge	RH		LH / RH
Door Handle	Lockable		
Type of Controller	Digital / Microprocessor		
Temperature Display	External		
Data Logger	Optional		
Door Ajar Alarms	Optional		
High and Low Readout	Visual / Audible		
BMS Connectivity			
(BMS Connectivity for Remote Alarms only)	Optional		
(BMS Connectivity for 4-20mA Transmitter)	Optional		
Port Entry	Optional		
Battery Back-up for Controller	N/A		
Power Requirement	10amp GPO		
Castors	Yes		
Adjustable Shelves	5		
Usable Shelf Area in mm(w/d)	480 x 450	750 x 530	1080 x 530
Light	LED		

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

# NLMi REFRIGERATOR INCUBATOR

**NULINE**  
REFRIGERATION

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to  $-40^{\circ}\text{C}$ ), blanket warmer, mortuary refrigerators and refrigerated incubators.

