
BLANKET WARMERS

NULINE REFRIGERATION

FEATURES

- Fan Forced System
- Overheating Alarm
- Safety Cutoff at 70°C
- Double Glazed Doors
- Stainless int./ext. with Chrome Plated Shelves
- Microprocessor Control
- Digital Display with High Alarm
- TGA Certified





700w x 700d x 920h 750w x 700d x 1880h 1200w x 700d x 1880h



BLANKET WARMERS



| Model | BW100 | BW260 | BW650 | BW1100 |
|--|-----------------|------------------|------------------|-------------------|
| Gross Capacity | 80 L | 247 L | 611 L | 1000 L |
| Blanket Capacity | 5 - 6 | 10 - 12 | 30 - 35 | 47 - 52 |
| External Dimensions in mm (w/d/h including castors) | 500 x 600 x 625 | 700 x 700 x 920 | 750 x 700 x 1880 | 1200 x 700 x 1880 |
| Internal Dimensions in mm (w/d/h) | 390 x 480 x 370 | 590 x 685 x 600 | 630 x 610 x 1590 | 1100 x 550 x 1555 |
| Packed Dimensions in mm (w/d/h including pallet) | 575 x 670 x 765 | 800 x 800 x 1070 | 850 x 800 x 2030 | 1300 x 800 x 2030 |
| Gross Weight - kg (packed weight) | 63 | 109 | 180 | 260 |
| Net Weight - kg | 48 | 89 | 160 | 220 |
| Voltage | 240V | | | |
| Temperature Range | 60°C | | | |
| Exterior Cabinet Finish | Stainless Steel | | | |
| Inner Liner Material | Stainless Steel | | | |
| Door | Glass | | | |
| Door Hinge | RH LH / RH | | | |
| Door Handle | Full Length | | | |
| Type of Controller | Microprocessor | | | |
| Temperature Display | Digital | | | |
| BMS Connectivity | | | | |
| (BMS Connectivity for Remote Alarms only) | Optional | | | |
| (BMS Connectivity for 4-20mA Transmitter) | Optional | | | |
| Port Entry | Optional | | | |
| Power Requirement | 10amp GPO | | | |
| Castors | Fixed Legs Yes | | | |
| Adjustable Shelves | 2 | 3 | | 6 |



Nuline Refrigeration reserves the right to alter the specifications without prior notice.





-

BLANKET WARMERS

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to -40°C), blanket warmer, mortuary refrigerators and refrigerated incubators.

