




ECF 289 2T COMBINATION

NULINE
REFRIGERATION

FEATURES

- Bottom mounted system
- Heated cabinet facia to prevent condensation
- Auto defrost in 2°C section
- Manual defrost in -40°C section
- Digital display with built in alarms
- Hi / Lo read out



	Refrigerator: 2°C to 8°C Freezer: -10°C to -40°C
	189 L – Refrigerator 100 L – Freezer
	700w x 640d x 1845h



ECF 289 2T COMBINATION



Model	Refrigerator	Freezer
Gross Capacity	189 L	100 L
External Dimensions in mm (w/d/h including castors)	700 x 640 x 1845	
Internal Dimensions in mm (w/d/h)	600 x 510 x 710	500 x 460 x 505
Packed Dimensions in mm (w/d/h including pallet)	830 x 750 x 2004	
Gross Weight - kg (packed weight)	165	
Net Weight - kg	144	
Voltage	220V	
Energy Consumption	2.5kWh/24h / @25°C	
Temperature Range	2°C to 8°C	10°C to -40°C
Ambient Temperature Range	16°C to 32°C	
Refrigerant	R600a	R290
Cooling System / Defrost	Fan Forced / Automatic	Static / Manual
Exterior Cabinet Finish	Enamelled Steel	
Inner Liner Material	Stainless Steel	
Door	Glass	Solid
Door Hinge	RH	RH
Door Handle	Pull	
Type of Controller	Digital	
Temperature Display	External	
Data Logger	Built-in	
Door Ajar Alarms	Yes	
High and Low Readout	Yes	
BMS Connectivity		
(BMS Connectivity for Remote Alarms only)	Yes	
(BMS Connectivity for 4-20mA Transmitter)	Optional	
Port Entry	Optional	
Battery Back-up for Controller	Yes	
Power Requirement	10amp GPO	
Watts	330W	
Castors	Yes	
Adjustable Shelves	3 Shelves + 1 drawer	3 drawers
Usable Shelf Area in mm(w/d)	1 x 585 x 370 2 x 585 x 418 drawer x 570 x 434 x 121	494 x 419 x 165 x
Light	LED	No