

Premium Vacc-Safe

VACCINE REFRIGERATORS



×



+



60 MONTH WARRANTY*

Premium Vacc-Safe

✉ sales@eurochill.com.au

🌐 www.eurochill.com.au

☎ 1300 459 140



VACC-SAFE 400 PREMIUM

VACC-SAFE 600 PREMIUM

VACCSAFE 1000 PREMIUM

- QCPP Compliant
- Independently Tested By The Pharmacy Guild Of NSW
- Operating Range 2-8°C
- Audible & Visual Hi/Lo Temp. Alarms
- Door Open Alarm
- Self-Closing Lockable Glass Door
- Touch Screen Customised Digital Controller
- High / Low Memory & Easy Reset Function
- White Powder Coated Steel Internal & External
- 4 Castors, with Two Front Locking
- External Access Port
- 24 Months Warranty*
- BMS External Alarm Contacts
- Bonus Medi Log 2 Product Simulation Data Logger Included With Audible & Visual Alarms

MODEL	VACC-SAFE 400 PREMIUM	VACC-SAFE 600 PREMIUM	VACC-SAFE 1000 PREMIUM
Capacity Litre	400	600	1000
Ext Dimensions Mm	606W * 620D * 1950H	700W * 790D * 2058H	1130W * 720D * 2023H
Int Dimensions Mm	513W * 495D * 1501H	603W * 665D * 1596H	1033W * 595D * 1561H
Shelves	7	6	12
Power Consumption	4.935kWh	6.209kWh	9.903kWh